




MINIMALLY INVASIVE SURGERY FOR P.F.C.[®] SIGMA[™] AND P.F.C.[®] SIGMA[™] RP KNEE SYSTEMS



SURGICAL TECHNIQUE

INTELLIGENT  ORTHOPAEDICS™

Contents

Ci™ System Connections	2	Tibial Resection Verification	35
Ci™ System Set-up	4	Registration of Epicondyles	36
Start-up Procedure	5	Registration of the Anterior Cortex	37
Ci™ System Menu Buttons	6	Whiteside's Line	38
Ci™ System Procedure Set-up	8	Femoral Condyle Modelling	39
Exposure	10	Anterior Cortex Modelling	40
Tibial and Femoral Array Positioning	11	Femoral Model Calculation and Verification	41
Tibial Array Fixation	12	Definition of the Anatomic Axis	42
Femoral Array Fixation	14	Alignment Settings	43
Camera Alignment	15	Ligament Balancing	44
Patellar Resection	16	Femoral Implant Planning	46
Registration Methods	18	Distal Femoral Resection Navigation	49
Ci™ Pointer Tips	20	Distal Femoral Resection Verification	51
Calculating the Femoral Head Centre	21	Anterior Femoral Resection Navigation	53
Femoral Mechanical Axis	22	Anterior Femoral Resection Verification	54
Definition of the Malleoli	23	Final Femoral Preparation	56
Definition of the Tibial Mechanical Axis	24	Tibial Preparation	59
Tibial Sizing	25	M.B.T. Tray Preparation	60
Tibial A/P Direction	26	M.B.T. DuoFix™ Tray Preparation	61
Tibial Modelling	27	Trial Reduction	62
Tibial Plateau Modelling	28	Final Component Implantation	64
Tibial Model Calculation and Verification	29	Final Patellar Implantation	65
Pre-operative Leg Alignment	30	Final Flexion and Extension Checks	67
Tibial Implant Planning	31	Instruments and Ordering Information	68
Tibial Resection	33		

Ci™ System Connections



Ci™ Unit Connections

Connecting the Ci™ System

All cables are connected to a panel at the rear of the Ci™ unit. All power cables and connectors are supplied with the system.

1. Power switch (on / off)

The power switch is positioned below the power socket. A green light indicates when the unit is switched on.

2. Power lead connection

The standard 3-pin power lead connects to the power socket.

3. Voltage switch

The red 'H-shape' switch adjusts the power supply between 110 volts and 220 volts. Make sure that the correct voltage is selected before switching on the system.

Power lead



Earth / ground connection

Camera cable



Camera unit cable 1



Camera unit cable 2



4. Earth / ground lead connection

The earth / ground socket is identified by the yellow and green colours.

The system must be earthed / grounded using the cable supplied before switching on the system.

5. Camera unit connection (tracking unit)

The orange coloured camera cable plugs into the camera unit socket (also marked in orange). The cable pins are aligned before inserting the cable. The cable is removed by pulling smoothly and slowly to avoid damaging the pins.


Important: The camera must be switched on at least 20 minutes prior to the start of the procedure.

Connections on the camera unit should also be checked at this stage.


Please refer to the System User Guide (Art-No: 52543-03EN) for detailed set-up instructions.

Note: The Ci™ instrumentation should be coated with a medical grade lubricant prior to surgery.

Ci™ System Set-up



Option 1: Ci™ Camera unit at the base of the table



Option 2: Ci™ Camera unit angled

Positioning the Ci™ System in the O.R.

The Ci™ unit should be positioned so that the surgeon has a clear view of the screen. The touch-screen should be within easy reach of the surgeon or surgical assistant.

Ideally the infrared camera unit should have an uninterrupted view of the tibial and femoral marker arrays throughout a full range of motion (knee flexion, extension and pivoting of the hip joint). The camera unit should be positioned between 1.5 to 2 metres (4 to 6 feet) away from the surgical site.

The Ci™ unit must be covered using the sterile drape supplied.

Unstable stands or tables should not be used for placement of the Ci™ unit.


Please refer to the System User Guide (Art-No: 52543-07EN) for detailed set-up instructions.

Start-up Procedure




Start-up screen

Figure 1



Camera connection error screen

Figure 2



System settings for network, video, language and touchscreen

Figure 3

All cables must be fully connected before switching on the system.

Once the power supply has been switched on, the Ci™ software will automatically initialise (Figure 1).


Two 'beeps' indicate that the camera unit is properly connected.


If the camera unit is not properly connected, no 'beeps' will be heard and an error message will appear asking you to restart the system (Figure 2).

The start-up screens allow the user to configure network and video settings if required (Figure 3).

A blank CD-R (not CD-RW) should be inserted into the Ci™ unit before draping to allow any patient data to be stored.

Ci™ System Menu Buttons





Ci™ System Procedure Set-up


On screen there are next and back buttons.
In this surgical technique we use arrows to indicate movement to the next step.

Arrow Key:
▶ Press 'next' to continue
◀ Press 'back' to return to the previous step



Patient name and ID input screen

Figure 4



Ci™ 1.1 or Ci™ MITKR workflow selection

Figure 5



Implant product line selection

Figure 6

The procedure is started by entering the patient's name and ID, using the touch-screen keyboard (Figure 4). This information will appear on the screen and on each screenshot taken during the operation. ▶

A choice of Ci™ Knee 1.1 or Ci™ MITKR software can be selected (Figure 5). ◀▶


The appropriate MITKR Knee System software is selected from the choices displayed - LCS® Complete™ / P.F.C.® Sigma™ / P.F.C.® Sigma™ RP Knee Systems (Figure 6). ◀▶



Navigation order selection

Figure 7

If the P.F.C.® Sigma™ product line is selected, a choice of 'tibia first' or 'optimised' navigation order options are presented (Figure 7). ◀▶



Specific implant selection

Figure 8

The appropriate product line is selected (Figure 8). Press 'change' to alter the selections on screen. ◀▶

The femoral component is selected first, followed by the tibial component and the tibial insert.



Affected knee selection

Figure 9

The appropriate treatment side for the surgery is selected (Figure 9). ◀▶

Exposure




Figure 10


Mid-vastus, sub-vastus or medial peripatellar capsular exposures are performed (Figure 10).

Following exposure of the joint, the anterior fat pad is cleared from the anterior cortex. The patella pouch and the working area around the tibia and femur should also be cleared.

In order to achieve long leg alignment and joint stability, it is vital that all osteophytes are removed from the entire medial or lateral edges of the femur and tibia before registration takes place.

A retractor is used on the patella rim to draw the patella into a central position. Any significant osteophytes on the edge of the patella should also be excised.

Tibial and Femoral Array Positioning



The arrays must be placed away from the intended positions of both the femoral and tibial instruments to avoid any contact during bone resection and drilling.

Note: The position of the reference arrays must not be moved during the operation as this will lead to inaccurate information being displayed.

Positioning of the cutting jig

The Ci™ MI cutting jig is held in an approximate position on the tibia before marking the pin positions. Care should be taken to position the array clear of the jig and any navigated instruments (Figure 11).

Positioning of the femoral array

The femoral array should be positioned away from the intended incision and should not influence soft tissue and ligament movement in extension and flexion (Figure 12).

Tibial Array Fixation



Figure 13

A fixation pin template is positioned on the tibia. Two marks are made through the template to identify the pin positions (Figure 13).




Figure 14

A stab incision is made through the soft tissue using a scalpel (Figure 14).




Figure 15

The two pins are fixed unicortically into the tibia. A sheath is used on the pins to avoid any disruption to the soft tissue. The array block is placed on the pins and secured in place with the locking screw (Figure 15).



Figure 16

The quick release array is secured to the block. Final positioning and tightening of the array is made using the adjustment screw (Figure 16).

Femoral Array Fixation



Figure 17

The Schantz screw position is marked on the femur and a stab incision is made through the soft tissue using a scalpel (Figure 17).



Figure 18

A single Schantz screw is drilled into the femur bicortically. A sheath is used on the screw to avoid any disruption to the soft tissue. The array sleeve is placed on the screw and pushed down until the spikes engage cortical bone.



The sleeve is secured using the locking and adjustment screws and tightened using the wrench (Figure 18).



Figure 19

The quick release array is secured to the block. Final positioning and tightening of the array is made using the array adjustment screw (Figure 19).

Camera Alignment



- The status field lights indicate that all instruments are visible
- Green: Normal tracking mode
 - Yellow: Ci™ instrument not visible to camera
 - Red: Ci™ instrument and / or navigation array not visible to camera for more than 3 seconds
 - Grey: Camera communication error during start-up
Re-start the system

Camera set-up screen

Figure 20

The camera alignment window highlights movement of the reference arrays and instruments. These are displayed as different coloured dots and circles on-screen (Figure 20).

Coloured dots represent the positions of the femoral and tibial reference arrays, and the Ci™ pointer.

The circles show the position of the instruments in relation to the camera unit's field of view.

The status field lights indicate the visibility of the arrays and instruments.

The camera displays can be accessed at any time during the procedure by touching the status field lights. ◀▶

Patellar Resection




Figure 21

The patella is everted to an angle of approximately 40 - 60 degrees. The synovial tissue and retinaculum are released from the periphery of the patella down to the plane of the quadriceps tendon.



Patella stylus

The thickness of the patella is measured and the level of bone resection calculated (Figure 21). The thickness of the resurfaced patella should be the same as the natural patella. The minimum residual dimension should be no less than 12 mm.




Figure 23

A patella stylus is selected that matches the depth of the patella to be resected (Figure 22). The leg is placed in extension with the patella cutting guide positioned with the stylus against the posterior cortex

of the patella with the serrated jaws at the superior and inferior margins of the articular surface (Figure 23).



Figure 24

The jaws should be closed to firmly engage the patella (Figure 24).




Figure 25

Figure 26

The stylus is removed and the resection performed using an oscillating saw with a 1.47 mm blade, through the saw capture and flush to the cutting surface (Figure 25).

A patella wafer is hand placed on the resected surface (Figure 26).

Registration Methods





Figure 27




Figure 28

Pivoting

Pivoting involves rotating the femur in a loose arc until a sufficient number of points have been acquired. The software uses this information to determine the rotational centre of the femoral head which is the most proximal point of the mechanical axis (Figure 27).

Single landmark acquisition

Key single landmarks are registered by touching the appropriate bone structure with the tip of the Ci™ pointer. On the monitor, a single arrow indicates where the point is to be acquired (Figure 28).

The following points are acquired in this way:

- Malleolus (medial and lateral)
- Mechanical axes (proximal tibia and distal femur)
- Proximal tibial contour (medial, lateral, anterior and posterior)
- Epicondyles (medial and lateral)
- Femoral anterior cortex resection level




Figure 29

Acquiring directions

Some axes and directions are acquired by holding the Ci™ pointer absolutely still in a specific alignment (Figure 29). The following directions are acquired in this way:

- Tibial rotation
- Whiteside's line

Multiple landmark acquisition

Multiple landmark acquisition registers bone areas. The points acquired are used to calculate resection levels and morph the 3D bone model. Points are acquired by sliding the Ci™ pointer tip along the bone structure (Figure 30).




Figure 30

Points are only acquired which are at least 1.5 mm from the previous five points.

Dots indicating the registered areas appear on the image of the bone with progress shown on the progress bar. The following points are acquired in this way:

- Tibial plateau (medial and lateral)
- Condyles (medial and lateral)
- Anterior cortex (tibia and femur)

Ci™ Pointer Tips




Figure 31



Figure 32

1. To acquire points accurately, it is helpful to hold the tip of the Ci™ pointer with one hand while pivoting the Ci™ pointer with the other. This will ensure that all of the points are registered on the bone, not away from the patient (Figure 31).
2. To avoid acquiring unnecessary points before the Ci™ pointer is positioned, it may be helpful to cover one of the marker spheres with a hand. This will effectively 'eliminate' the Ci™ pointer from the camera's field of view (Figure 32).
3. Make sure all three marker spheres on the Ci™ pointer are directed towards the camera unit when acquiring a point.
4. When a point has been successfully acquired, the system will sound an audible 'beep'. This removes the need to refer to the screen as each point is acquired. If the camera cannot see the Ci™ pointer, a deep sounding 'beep' will be heard.
5. If at any stage during the procedure the marker spheres become contaminated with blood or tissue, they may be gently cleaned using lint-free cloth.

Calculating the Femoral Head Centre



An assistant may hold the pelvis in position during this step to ensure that an accurate calculation is made. Sidebars on the OR table may inhibit the registration process.

Figure 33



Figure 34

Calculating the rotational centre of the femoral head defines the proximal point of the mechanical axis. It is important to make sure that the patient's pelvis is not moved excessively during registration as this will lead to miscalculation by the system.

The femur is pivoted using circular motions (Figure 33). The system calculates a series of points to determine the rotational centre and will automatically proceed when the rotational centre has been accurately calculated (Figure 34). ◀▶

Pivoting will restart automatically if the error is more than 3 mm.

After the fourth attempt (and if the error is between 3 mm and 5 mm) the error value is displayed on screen. The surgeon has the option to accept the error or try again.

If an error of 4 mm or 5 mm is accepted, the mechanical axis will deviate by more than 0.5 degrees.

Femoral Mechanical Axis



Figure 35

The Ci™ pointer should be placed slightly medial to the anterior aspect of the femoral notch point (as indicated on screen) (Figures 35 & 36).

Do not use the femoral canal entry point as a guide for acquiring the point. The acquisition of this point completes the femoral mechanical axis. ◀▶




Figure 36

Definition of the Malleoli



Figure 37





Figure 38



The medial and lateral malleoli are defined using the Ci™ pointer.

The malleoli can usually be located by hand before acquiring the points. It is important that draping or bandaging is reduced to a minimum to enable the malleoli to be located.

Place the tip of the Ci™ pointer on the medial malleolus and pivot the Ci™ pointer (Figure 37). It is important not to move the tip of the Ci™ pointer from the bone. Once the system has registered the medial malleolus, the lateral malleolus can be registered in the same way (Figure 38).

Acquiring the malleoli defines the distal point of the mechanical axis.

All single registration points for anatomic landmarks are acquired in the same way. ◀▶

It is important that these points are acquired as accurately as possible.

Definition of the Tibial Mechanical Axis



Figure 39

The final point on the mechanical axis is defined by acquiring the posterior aspect of the ACL tibial insertion point (Figure 39). This is indicated by the arrow on screen (Figure 40).

It is important that the point acquired defines the mechanical axis, and not the tibial eminence.




Figure 40

The definition of the mechanical axis is the basis for all further calculations and should be acquired as accurately as possible. Final implant position will be referenced to the mechanical axis. ◀▶

Tibial Sizing



Figure 41

After removal of any osteophytes, the most medial and lateral points of the proximal tibia are acquired using the Ci™ pointer (for implant sizing). These are followed by the most anterior and posterior points (for implant positioning) (Figures 41 & 42).

The points should be registered at the expected resection level. This will help to avoid oversizing the tibial component. ◀▶







Figure 42

Note: Failure to remove osteophytes at this stage may lead to oversizing the implant (during implant planning).

Tibial A/P Direction



Figure 43

The Ci™ pointer is placed on the tibia to acquire the A/P direction only (Figure 43). If access to the centre of the tibia is limited, the Ci™ pointer can be placed to one side. Accurate acquisition of the A/P direction will help to avoid a compound slope.

The Ci™ pointer is held in place for 2 seconds to allow the system to calculate the A/P direction (Figure 44). ◀▶




Figure 44

Tibial Modelling




Figure 45

A fixed number of points along the anterior and medial / lateral tibia are acquired using the Ci™ pointer. These points are used to further define the bone model (Figure 45).

The Ci™ pointer is placed on the tibia and pivoted to begin the procedure. A 'beep' will indicate when to begin moving the Ci™ pointer (Figure 46).




Figure 46

The tip of the Ci™ pointer should 'paint' the exposed surface of the tibia making sure it does not leave the surface of the bone. Points should be acquired mainly around the rim and extend as posteriorly as possible. Points further down the tibia

are less important, but will add shape to the model. The system will indicate when enough points have been acquired. ◀▶

Tibial Plateau Modelling



Figure 47

The Ci™ pointer should be placed in the centre of the medial tibial plateau. Laminar spreaders can be used to open the joint to allow greater access to the tibial plateau. The tip of the Ci™ pointer is moved in a circular motion across the surface of the plateau (Figure 47).





Figure 48

It may be helpful to apply varus or valgus stress to the knee joint to allow greater access to the tibial plateau.

The system will automatically proceed to the lateral plateau, which is registered in an identical manner to that of the medial plateau (Figure 48). ◀▶

Tibial Model Calculation and Verification




Figure 49

Following registration, the system adapts the tibial model (Figure 49).

- White areas are parts of the model with deviations no greater than 2 mm from the points acquired during registration.

- Brown areas are parts of the model where no points have been acquired.

The accuracy of the model is checked by holding the Ci™ pointer to the tibia. The exact deviation from the tip of the Ci™ pointer to the model is displayed




Figure 50

on the screen. The maximum acceptable deviation is less than 1 mm. Acquired points will normally show a deviation of less than 1 mm. Verification of the model will only be accurate in white areas where points have been acquired (Figure 50).

Note: Implant sizing is based on acquired points only. The accuracy of the model should be checked on the areas responsible for implant sizing and position (medial and lateral boundary, medial and lateral plateau).

Pre-operative Leg Alignment




Figure 51

It is important to register the knee joint in extension. The leg is held in clinical extension for 2 seconds to register the patient's pre-operative condition prior to the procedure (Figure 51).

- The orange lines on screen represent the tibial and femoral mechanical axes.
- The green line represents the hip / knee / ankle (HKA) line.





Figure 52

This step allows a comparison to be made between the initial and final alignments at the end of the procedure (Figure 52). ◀▶


Tibial Implant Planning



Tibial implant planning screen

Figure 53

Initial planning and alignment of the implant position is automatically calculated during the registration procedure. The system will select the thinnest implant available based on the points acquired.




Varus / valgus rotation

Resection level

Varus / valgus rotation

Figure 54

Fine tune functions allow the user to reposition the implant as required. Tibial slope, varus / valgus and resection level can be adjusted on screen at this stage (Figures 53, 54 & 55).



Slope adjustment

Slope adjustment

Figure 55

Button Functions

Blue arrows – allow repositioning of the tibial component.

Yellow line – tibial resection level.

The image can be enlarged on screen using the zoom function on the menu bar.




Figure 56

The A/P position of the implant is aligned to the tibial anterior point determined during registration (Figure 56).

The following information is displayed in the bottom right hand side of the screen (Figure 57):

Posterior slope – (in degrees).

Resection high - indicates the distance of the deepest point on the higher plateau from the resection level (in mm).




Figure 57

The + and - buttons change the size of the tibial component. The tibial cut is re-calculated and the gaps adjusted accordingly. The system only allows compatible sized implants to be selected.

Resection low - indicates the distance of the deepest point on the lower plateau from the resection level (in mm).

Varus / valgus – (in degrees). ◀▶

Tibial Resection




Figure 58

Before placing the Ci™ MI cutting jig on the tibia, its varus / valgus, A/P and proximal translation positions should be set as close to a mid-position as possible to allow for fine-tune adjustment in all directions. The grooves on the adjuster threads should be used as a reference for the mid-position.

The insertion plate on the navigated Ci™ cutting block adaptor is placed into the saw capture on the Ci™ MI cutting jig (Figure 58). The yellow plane represents the planned plane and the blue plane represents the current plane (Figure 59).




Figure 59

The screen data is used to guide the jig into an approximate position (in line with the mechanical axis of the tibia and close to the desired tibial resection level).

The jig is fixed in place on the tibia using a Schanz screw. The locking trigger is pulled forward to secure the jig to the screw. ◀▶




Figure 60

The varus / valgus alignment of the instrument should be checked and adjusted if necessary. When this is complete, a second pin is positioned in the middle of the slot to support the instrument in this position (Figure 60).

The tibial slope, varus / valgus and proximal translation adjustments can now be made using the three adjustment screws, with reference to the on screen data.



Figure 61

Once accurately positioned, the jig is secured using extra pins below the saw capture (Figure 61).

The navigated Ci™ cutting block adaptor is removed and the tibial resection completed with a 1.47 mm saw blade.

The tibial resection must always be made through the saw capture slot.

The leg is placed in approximately 10-15 degrees of flexion and the resected tibia is removed.

Tibial Resection Verification



Figure 62

The insertion plate of the Ci™ cutting block adaptor is placed flat on the resected plane (Figure 62).

Any difference between planned and current resection planes can be seen on screen (Figure 63).




Figure 63

The yellow plane represents the planned plane and the blue plane represents the current plane. ◀▶




Figure 64

The tibial resection numeric data is displayed by pressing 'next'. The tibial resection can be repeated if necessary (Figure 64). ◀▶

Registration of Epicondyles




Figure 65





Figure 66

Registration of the femur begins once the tibial resection has been verified.

Acquisition of the most medial and lateral epicondylar points are used to define the epicondylar axis (Figures 65 & 66). ◀▶

Registration of the Anterior Cortex



Figure 67

The anterior cortex is acquired using the Ci™ pointer. Ideally, the Ci™ pointer should be placed on a point on the anterior femoral cortex just above the superior border of the patello-femoral articular surface (Figure 67).

It should indicate the place where the anterior flange of the femoral component would sit if flush to the anterior cortex.




Figure 68

This will enable the system to determine the size of the femoral implant required (Figure 68). ◀▶

It is important to acquire this accurately to avoid the intended cutting plane notching the femur.

Whiteside's Line




Figure 69



Figure 70




Figure 71

Whiteside's line is used as an optional reference for femoral component rotation alignment. It should be marked initially using electro-cautery. It is easiest to draw by looking along the horizon of the trochlear groove. Once the line is drawn the Ci™ pointer can be held along this line.

The Ci™ pointer must be held perfectly still while the system acquires this reference (Figures 69, 70 & 71). In the event of a system failure the marked line will be used as an alignment reference point for the back-up instrument set. ◀▶

Femoral Condyle Modelling




Figure 72

A fixed number of points along the surface of the medial and lateral condyles are acquired using the Ci™ pointer (Figures 72 & 73).





Figure 73

The tip of the Ci™ pointer should 'paint' the surface of the condyles. Points should be acquired as posteriorly as possible and along the distal part of the affected condyles. This allows the Ci™ System to accurately calculate both the size and distal resection level of the intended implant.

If insufficient posterior points are acquired the system will ask for the step to be repeated. ◀▶

Anterior Cortex Modelling




Figure 74




Figure 75

Multiple points along the anterior cortex are acquired using the Ci™ pointer. These points are used to define the surface and to help determine if the implant is likely to notch the femur (Figures 74 & 75). ◀▶

Femoral Model Calculation and Verification




Figure 76

Following registration, the system adapts the femoral model (Figure 76).

- White areas are parts of the model with deviations no greater than 2 mm from the points acquired during registration.

- Brown areas are parts of the model where no points have been acquired.

The accuracy of the model is checked by holding the Ci™ pointer to the tibia. The exact deviation from the tip of the Ci™ pointer to the model is displayed




Figure 77

on the screen. The maximum acceptable deviation is less than 1 mm. Acquired points will normally show a deviation of less than 1 mm. Verification of the model will only be accurate in white areas where points have been acquired (Figure 77).

Note: Femoral implant sizing is based on acquired points only. The accuracy of the model should be checked on the areas responsible for implant sizing and position.

Definition of the Anatomic Axis (Femoral Bow)





Figure 78


The flexion / extension and rotation position of the femoral implant can be positioned according to the anatomic or mechanical axis. If the mechanical axis is required, continue to the next step. If the anatomic axis is required the Ci™ cutting block adaptor is placed in the yoke on



Anatomic axis (femoral bow) acquisition screen

Figure 79

the anterior femur and held still for 2 - 3 seconds to register the anatomic axis (Figure 78). A blue plane is displayed on the screen indicating the position of the Ci™ cutting block adaptor (Figure 79). ◀▶




Anatomic and mechanical axis positions

Figure 80

The software will display the angular difference between the mechanical and anatomic axes in the sagittal plane (Figure 80). The femoral component position will be rotated by the same angle to ensure accurate anterior alignment. ◀▶

If a negative angle is recorded the implant will be positioned according to the mechanical axis to avoid femoral notching.

Alignment Settings (P.F.C.® Sigma™ Knee System Tibia First - Non-Optimised Workflow Only)



Implant alignment and rotation reference settings

Figure 81




Figure 82

If the tibia first workflow was selected at the start of the procedure, a choice of anterior or posterior alignment options are displayed for positioning the femoral implant. Rotation can be referenced to either the epicondylar line, posterior line, or Whiteside's line (Figures 81 & 82).

Ligament Balancing





Figure 83

Ligament balancing allows the ligament tension on the medial and lateral collateral ligaments to be assessed. A tensioning device such as laminar spreaders, knee tensioners or Moreland knee retractors should be used to balance the joint.

Ligament balancing in extension
The leg is placed in extension and distracted. Varus and valgus stresses are applied to the knee joint. The laminar spreaders can be used to achieve this (Figure 83).

Ligament balancing in extension



Information button

Figure 84

Pressing the information button displays the gap information and additional angular values

The Ci™ System displays the varus / valgus and flexion angle. Pressing the information button will display the gap information and additional angular values.

The values for the extension gap, varus / valgus angle and flexion angle are stored on the system by pressing the 'store' button and are used as the basis for femoral planning (Figure 84).




Figure 85

Ligament balancing in flexion

The leg is placed in flexion. Independent but equal pressure is applied to the medial and lateral compartments of the knee (Figure 85).

Pressing the information button will display the gap information and additional angular values.

Ligament balancing in flexion



Information button


Figure 86

Pressing the information button displays the gap information and additional angular values

The values for the flexion gap, varus / valgus angle and flexion angle are stored on the system by pressing the 'store' button (Figure 86).

In the optimised workflow this data will be used to rotate the implant with respect to the proximal tibial resection in order to achieve a parallel gap space in flexion. ◀▶

Femoral Implant Planning



First femoral implant planning screen

Figure 87

The data gathered during the registration and verification procedures is used to calculate the femoral implant size, position and resection levels. If the optimised workflow was selected, the system will size and position the femoral component in order to achieve equal and parallel flexion / extension gaps. If the tibia first workflow

was selected the flexion / extension position may need to be manually adjusted. The gap values will only be available if they were stored during ligament balancing. Any recommendations made by the system should be thoroughly checked and adjusted as necessary.

The rotation of the femoral implant is based around the anterior cortex point selected earlier (shown by red dot on anterior cut line). The A/P shift data is defined using this point as a reference




Figure 88

Button functions and screen information:

Blue arrows – allow repositioning of the femoral component.

Dark blue lines – proposed cutting planes of the femur.

Orange line - mechanical axis line.

The points acquired during registration provide the true reference for each of the fine tune functions. The 3D model only provides additional information and 3D orientation (Figure 87). The following adjustments to the plan can be made at this stage: flexion / extension, resection level and A/P shift (Figure 88).

Planning information



The implant size / freesize function is used to change the size of the femoral component.

Resets all adjustments

Returns to the previous step

Accepts all adjustments

Figure 89

If the extension / flexion gap is too small for the thinnest available insert, the measurement(s) will be displayed in red. The minimum amount of extra bone resection will be displayed in millimetres (Figure 89).

The '+' and '-' buttons are used to change the femoral implant size. The 'freesize' option can be used to manually change both the tibial and femoral implant size. However, if inappropriate implants are selected, the system will display a warning message relating to implant compatibility.

Freesize warning screen





Figure 90

The 'proceed' option will be removed at this stage to force the user to select compatible implants (Figure 90).



Second femoral implant planning screen

Figure 91

The second femoral planning screen allows any of the following adjustments to the femoral component to be made if required: M/L position, varus / valgus and internal / external rotation (Figure 91).




Figure 92

The M/L position can be adjusted on screen with the blue arrows (a neutral position is recommended by the system – central to the mechanical axis).

Varus / valgus can be adjusted with the blue arrows. Rotational adjustment is made around the most distal point on the condyle (indicated by a blue dot).

The internal / external rotation value is set in relation to the chosen alignment (epicondylar / posterior / Whiteside's line). Deviation from other rotational references are displayed on screen (Figure 92).

Distal Femoral Resection Navigation




Figure 93

The Ci™ MI cutting jig is rotated 180° and positioned on the femur. The varus / valgus, flexion / extension and internal / external rotation positions should be set as close to neutral as possible. The grooves on the adjuster threads should be used as a reference for the mid-position.

The Ci™ cutting block adaptor is placed in the saw capture slot of the cutting jig and the assembly is approximately navigated into position on the femur using the on-screen data (Figures 93 & 94).

An initial pin is introduced to fix the jig and the lock is applied.




Figure 94

The varus / valgus alignment of the Ci™ MI cutting jig should be checked and adjusted if necessary. ◀▶




Figure 95

With the Ci™ MI cutting jig in the correct position, a second pin is placed in the middle of the slot to support the instrument.

Final fine tuning of varus / valgus, flexion / extension and distal femoral resection level can be carried out using the three adjustment screws (Figure 95). The Ci™ MI cutting jig is secured using additional pins in to the saw capture block.



Figure 96

The navigated Ci™ cutting block adaptor is removed and the distal femoral resection completed with the knee in extension using a 1.47 mm saw blade (Figure 96).

The distal femoral resection must always be made through the saw capture slot.

Distal Femoral Resection Verification




Figure 97

The insertion plate of the Ci™ cutting block adaptor is placed flat on the resected plane (Figure 97).





Figure 98

Any difference between the planned and current resection planes can be seen on screen (Figure 98). ◀▶




Distal femoral resection verification screen

Figure 99

The yellow plane represents the planned cut.

The blue plane represents the current cut.



Data showing the difference between the planned resection and the current resection

Figure 100

The resection numeric data is displayed by pressing 'next' (Figures 99 & 100). ◀▶

The resection can be repeated if necessary.

Anterior Femoral Resection Navigation




Figure 101

The appropriate side anterior resection cutting guide is selected (left or right) and positioned flush to the resected distal femur. The cutting guide should be set to the mid position (size 3), using the scale as a reference (Figure 101).




Figure 102

The Ci™ cutting block adaptor is placed into the slot on the cutting guide. The screen data is used to approximately position the instrument in the correct place.




Figure 103

The initial pin is introduced to fix the cutting guide on the resected femur. Rotation is navigated using the on-screen data (Figure 103). When complete, the fixation button is pressed to secure the lower half of the jig (Figure 102).

The resection height is navigated using the adjustment screw and the cutting guide is fixed in place with a final pin (Figures 102 & 103). ◀▶

Anterior Femoral Resection Verification




Figure 104

The anterior femoral resection is completed using an oscillating saw with a 1.47 mm blade (Figure 104).





Figure 105

The insertion plate of the Ci™ cutting block adaptor is placed flat on the resected surface (Figure 105).




Anterior femoral resection verification screen

Figure 106

Any difference between the planned and current resection planes can be seen on screen (Figure 106). ◀▶

The yellow plane represents the planned cut.
The blue plane represents the current cut.



Data showing the difference between the planned resection and the current resection

Figure 107

The resection numeric data is displayed by pressing 'next' (Figure 107). ◀▶

The resection can be repeated if necessary.

Final Femoral Preparation




Figure 108

The appropriate sized finishing block is selected and positioned on the resected anterior and distal surfaces of the femur.

The fixation lugs on each side of the block correspond exactly to the M/L dimension of the final implant and are used to visually place the instrument in the correct M/L position (Figure 108).




Figure 109

The block is pinned in place through the fixation lugs before any bone cuts are made.

Final resections are made using a 1.47 mm saw blade in the following order: Posterior femoral (Figure 109), chamfer cuts, box cut / peg holes.




Figure 110




Figure 111

The chamfer cuts are made with the saw flush to the angled plates on the cutting block (Figures 110 & 111).




Figure 112




Figure 113

The box cut is made using an oscillating or reciprocating saw. An osteotome may be used to clear any remaining bone (Figures 112 & 113).

Tibial Preparation




Figure 114

The tray trial allows for standard, M.B.T., keeled and 4 peg tray preparation (Figure 114). The MI reamer tower is attached to the tray trial by lifting the catch and rotating the tower until the button engages in the tray (Figure 115).

The tower is locked in place when the catch is released. The tibial reaming depth is controlled by a series of modular collars attached to the reamer tower.



Figure 115

With the correct collar attached, the tibia is reamed to the appropriate depth.



Figure 116

The stem punch is advanced into the tower and impacted into the cancellous bone until the appropriate tray size marking is reached (Figure 116).

M.B.T. Tray Preparation




Figure 117




Figure 118



Figure 119

Keeled Tray Option

If a keeled M.B.T. tray is to be employed, and the bone of the medial or lateral plateau is sclerotic, it is helpful to initially prepare the keel slot with an oscillating saw or high speed burr.

The universal handle is assembled to the appropriately sized M.B.T. keel punch and inserted into the M.B.T. punch bushing, taking care to avoid malrotation.

The assembly is impacted into the cancellous bone until the shoulder of the punch is in even contact with the M.B.T. punch bushing (Figures 117 & 118).

Non-Keeled Tray Option

The RP adaptor is inserted into the punch bushing (Figure 119).

M.B.T. DuoFix™ Tray Preparation



Figure 120

To finalise the rotational position, the four peg holes must be prepared. Remove the insert trial, femoral trial and APG or RP adaptor. Attach the tibial alignment guide handle, and verify that orientation of the tibial trial corresponds to the markings on the tibia (Figure 120).



Figure 121

Pin the tibial template, and drill the four peg holes (Figure 121).

Trial Reduction




Figure 122
(cruciate retaining)

The trial bearing size determined during implant planning is selected and inserted onto the M.B.T. tray trial.




Figure 123
(cruciate substituting)

The knee is placed in deep flexion. The femoral trial is inserted onto the distal femur (Figures 122 & 123).



Motion analysis with trial components in extension

Figure 124



Motion analysis with trial components in flexion

Figure 125

With the trial components in place, the tracking pattern is assessed through a full flexion cycle. The Ci™ pointer should be held along the tibial tray to check rotation. If the position of the Ci™ pointer matches the motion analysis pattern, the rotation position of the tray should be marked on the tibia (Figures 124 & 125).

The following points should also be assessed:

- Adequate range of motion
- Medial and lateral lift-off (using gap data)

- Proper ligamentous tension in extension and in flexion
- Correct mechanical alignment of the extremity
- Natural motion without restrictions

The green dashed line represents the registered A/P direction of the tibial bone. The solid red line represents the A/P direction of the femoral component. The red dots mark femoral/tibial contact through the flexion cycle. The pattern can be cleared by pressing 'reset'.

Final Component Implantation



Figure 126
(cruciate retaining)



Figure 127
(cruciate substituting)

The components are implanted in the following order (Figures 126 & 127).

Tibial Implantation

The entire site should be thoroughly cleansed using pulsatile lavage. Bone cement is prepared and applied by syringe or with digital pressure in its low viscous state to assure maximal penetration into the trabecular bone.

The universal handle is attached to the tray inserter. The tibial tray is carefully inserted, avoiding malrotation. When fully inserted, several mallet blows may be delivered to the top of the universal handle. All extruded cement is removed using a curette.

Polyethylene Implantation

Loose fragments or particulates should be removed from the permanent tibial

tray. The appropriate permanent tibial insert can be inserted at any time during the cementing procedure.

Femoral Implantation

The entire site should be thoroughly cleansed using pulsatile lavage. Bone cement is prepared and applied to the femur. The femur is hyperflexed and the tibia is subluxed forward. The implant

is attached to the femoral inserter and inserted onto the femur. The knee is extended to approximately 90° for final impaction. The inserter is released and any extruded cement is cleared using a curette.

Final Patella Implantation



Figure 128

A template is selected that most adequately covers the articular surface without overhang. The patella wafer is removed from the patella and the handle positioned on the medial side of the everted patella.

Where bone is deficient on the lateral side, the next smaller size is selected and positioned slightly medially to enhance patellar tracking (Figure 128).



Figure 129

The template is firmly engaged to the resected surface and the holes made with the appropriate drill bit (Figure 129).




Figure 130



Figure 131


The patellar implant may now be cemented. The cut surface is thoroughly cleansed with pulsatile lavage. Cement is applied to the surface and the component is inserted.

The patellar clamp is designed to fully seat and stabilise the implant as the cement polymerises. It is positioned with the silicon O-ring, centred over the articular surface of the implant and the metal backing plate against the anterior cortex, avoiding skin entrapment.

When snug, the handles are closed and held by the ratchet until polymerisation is complete. All extruded cement is removed with a curette.


The clamp is released by unlocking the LOCK switch and squeezing the handles together (Figure 130). The patella is reduced and the patella implant is evaluated. An unrestricted range of motion, free bearing movement and proper patellar tracking should be evident (Figure 131).

Final Extension and Flexion Checks




Motion analysis with final components in extension

Figure 132



Motion analysis with final components in flexion

Figure 133



Shut down screen

Figure 134

Final motion analysis checks are made in both extension and flexion to assess joint kinematics (Figures 132 & 133). ◀▶

All planning information, screenshots and therapy reports can now be recorded to the CD-R.

This information can be saved prior to shut-down. The system shut-down begins after pressing the 'close' button on the final screen.

On the close application screen, press 'finish' to exit the Ci™ software (Figure 134). ◀

Instruments and Ordering Information

MI 8200	Base Kit MI CAS	950501519	MITKR CAS M.B.T. Evaluation Bullet Size 4-7
950501543	MITKR CAS Tibial / Femoral Jig	950501523	MITKR CAS RP Tibial Bullet Size 2-3.5
950501545	MITKR CAS 0 Degree 55 mm Cutting Block	950501524	MITKR CAS RP Tibial Bullet Size 4-7
MI 8300	P.F.C.® Sigma™ Femoral MI CAS	966520	Universal Handle
950501541	MITKR CAS Anterior Jig Sigma™ Assembly	217830118	M.B.T. Central Drill
950501528	MITKR CAS Femoral 3-in-1 Sigma™ Size 2	217830119	M.B.T. Central Stem Punch
950501529	MITKR CAS Femoral 3-in-1 Sigma™ Size 2.5	250025000	M.B.T. Modular Depth Stop Ring
950501530	MITKR CAS Femoral 3-in-1 Sigma™ Size 3	217830121	M.B.T. Plateau Trial Post
950501531	MITKR CAS Femoral 3-in-1 Sigma™ Size 4	217830137	M.B.T. RP Trial Button
950501532	MITKR CAS Femoral 3-in-1 Sigma™ Size 5	965383	M.B.T. tray impactor
MI 8400	M.B.T. MI CAS	MI 8500	LCS® Femoral MI CAS
950501502	MITKR CAS Tibial Template Size 2	950501542	MITKR CAS Anterior Jig LCS® Assembly
950501503	MITKR CAS Tibial Template Size 2.5	950501535	MITKR CAS Femoral 4-in-1 LCS® Small+
950501504	MITKR CAS Tibial Template Size 3	950501536	MITKR CAS Femoral 4-in-1 LCS® Medium
950501506	MITKR CAS Tibial Template Size 4	950501537	MITKR CAS Femoral 4-in-1 LCS® Standard
950501508	MITKR CAS Tibial Template Size 5	950501538	MITKR CAS Femoral 4-in-1 LCS® Standard+
950501511	MITKR CAS Bayonet Tower Assembly	950501539	MITKR CAS Femoral 4-in-1 LCS® Large
950501512	MITKR CAS Bushing Extractor	950501547	MITKR CAS 15 Degree 55 mm Cutting Block
950501514	MITKR CAS KL Punch Bushing Size 2-3	MI 8600	P.F.C.® Sigma™ RP-F MI CAS
950501515	MITKR CAS KL Punch Bushing Size 4-7	950501552	MITKR CAS Femoral 3-in-1 RP-F Size 2
950501516	MITKR CAS M.B.T. Bullet WO Pegs Size 1-3	950501553	MITKR CAS Femoral 3-in-1 RP-F Size 2.5
950501517	MITKR CAS M.B.T. Bullet WO Pegs Size 4-7	950501554	MITKR CAS Femoral 3-in-1 RP-F Size 3
950501518	MITKR CAS M.B.T. Evaluation Bullet Size 1-3	950501555	MITKR CAS Femoral 3-in-1 RP-F Size 4

950501556 MITKR CAS Femoral 3-in-1 RP-F Size 5
 954120 RP-F Tibial Insert Trial 10 mm Size 2
 954121 RP-F Tibial Insert Trial 12.5 mm Size 2
 954122 RP-F Tibial Insert Trial 15 mm Size 2
 954123 RP-F Tibial Insert Trial 17.5 mm Size 2
 954125 RP-F Tibial Insert Trial 10 mm Size 2.5
 954126 RP-F Tibial Insert Trial 12.5 mm Size 2.5
 954127 RP-F Tibial Insert Trial 15 mm Size 2.5
 954128 RP-F Tibial Insert Trial 17.5 mm Size 2.5
 954130 RP-F Tibial Insert Trial 10 mm Size 3
 954131 RP-F Tibial Insert Trial 12.5 mm Size 3
 954132 RP-F Tibial Insert Trial 15 mm Size 3
 954133 RP-F Tibial Insert Trial 17.5 mm Size 3
 954140 RP-F Tibial Insert Trial 10 mm Size 4
 954141 RP-F Tibial Insert Trial 12.5 mm Size 4
 954142 RP-F Tibial Insert Trial 15 mm Size 4
 954143 RP-F Tibial Insert Trial 17.5 mm Size 4
 954150 RP-F Tibial Insert Trial 10 mm Size 5
 954151 RP-F Tibial Insert Trial 12.5 mm Size 5
 954152 RP-F Tibial Insert Trial 15 mm Size 5
 954153 RP-F Tibial Insert Trial 17.5 mm Size 5
 954212 RP-F Trial Femur Size 2 Left
 954213 RP-F Trial Femur Size 2.5 Left
 954214 RP-F Trial Femur Size 3 Left

954215 RP-F Trial Femur Size 4 Left
 954216 RP-F Trial Femur Size 5 Left
 954222 RP-F Trial Femur Size 2 Right
 954223 RP-F Trial Femur Size 2.5 Right
 954224 RP-F Trial Femur Size 3 Right
 954225 RP-F Trial Femur Size 4 Right
 954226 RP-F Trial Femur Size 5 Right

MI 7400

966554 Base Tibial Instruments
 966341 Tibial Tray Trials Size 2
 966343 Tibial Tray Trials Size 3
 966344 Tibial Tray Trials Size 4
 966345 Tibial Tray Trials Size 5
 860590 Trial Tray Fixation Pins (pack of 4)
 966430 Tibial Trial Insert Remover
 966520 Universal Handle
 966342 Tibial Tray Trial Size 2.5

SP 2714

P.F.C.® Sigma™ Fixed Bearing MI CAS
 Specialist® 2 Standard Tibial
 966350 Standard Tibial Punch Guide Size 1.54
 966351 Standard Tibial Punch Guide Size 5 & 6
 966355 Standard Non Cemented Tibial Punch Size 1.52
 966356 Standard Non Cemented Tibial Punch Size 2.54
 966357 Standard Non Cemented Tibial Punch Size 56

966358 Standard Cemented Tibial Punch Size 1.52
966359 Standard Cemented Tibial Punch Size 2.56
966564 Cruciform Keel Punch Sterile Tray Insert

SP 2715 Specialist® 2 Modular Cemented Punch Instruments

966360 Modular Tibial Cemented Punch Guide Size 1.5-3
966361 Modular Tibial Cemented Punch Guide Size 4-6
966370 Modular Tibial Cemented Punch Size 1.5-3
966371 Modular Tibial Cemented Punch Size 4-5
966376 Specialist® 2 Modular Tibial Drill Bush Cemented Size 1.5-3 (13 mm)
966378 Specialist® 2 Modular Tibial Cemented Drill Bush Size 4-5 (15 mm)
966380 Modular Tibial Drill Size 1.5-3 13 mm Cemented
966381 Modular Tibial Drill Size 4-6 15 mm Cemented
966566 Cemented Modular Punch

MI 8800 Demo MI CAS (International Only)

950501504 MITKR CAS Tibial Template Size 3
950501511 MITKR CAS Bayonet Tower Assembly
950501512 MITKR CAS Bush Extractor
950501514 MITKR CAS Keel Punch Bush Size 2-3
950501516 MITKR CAS M.B.T. Bullet WO Pegs Size 1-3
950501518 MITKR CAS M.B.T. EVAL Bullet Size 1-3
950501521 MITKR CAS DuoFix™ Tibial Reamer
950501523 MITKR CAS RP Tibial Bullet Size 2-3.5
950501541 MITKR CAS Anterior Jig Sigma™ Assembly

950501542 MITKR CAS Anterior Jig LCS® Assembly
950501530 MITKR CAS Femoral 3-in-1 Sigma™ Size 3
950501537 MITKR CAS Femoral 4-in-1 LCS® Standard
950501554 MITKR CAS Femoral 3-in-1 RP-F Size 3
950501543 MITKR CAS Tibial/Femoral Jig
950501545 MITKR CAS 0 Degree 55 mm Cutting Block
950501547 MITKR CAS 15 Degree 55 mm Cutting Block

MI 8000 LCS® Mobile Bearing MI CAS Instruments

180118000 1/8 Inch Drill Bit
228769000 Calliper
217856025 LCS® Completion Femoral Impactor/Extractor
864348 Bone File
966515 Specialist® 2 Pin Puller
966520 Universal Handle
869189 P.F.C.® J&J Tibial Retractor
217842005 Completion™ QDrill Driver
217842010 Completion™ QDrill Pin Pack
228748000 Femoral Lug Drill
258919000 Visualisation Wing
966180 Modular Femoral Impactor Head

Individual codes

950501500	MITKR CAS Tibial Template Size 1	950501577	MITKR CAS Patella Trial - LCS® Large
950501501	MITKR CAS Tibial Template Size 1.5	950501578	MITKR CAS Patella Trial - LCS® Large+
950501505	MITKR CAS Tibial Template Size 3.5	950501579	MITKR CAS Patella Trial - Sigma™ 25 mm
950501507	MITKR CAS Tibial Template Size 4.5	950501580	MITKR CAS Patella Trial - Sigma™ 28 mm
950501509	MITKR CAS Tibial Template Size 6	950501581	MITKR CAS Patella Trial - Sigma™ 32 mm
950501510	MITKR CAS Tibial Template Size 7	950501582	MITKR CAS Patella Trial - Sigma™ 35 mm
950501513	MITKR CAS Keel Punch Bush Size 1-1.5	950501583	MITKR CAS Patella Trial - Sigma™ 38 mm
950501522	MITKR CAS RP Tibial Bullet Size 1-1.5	950501584	MITKR CAS Patella Trial - Sigma™ 41 mm
950501527	MITKR CAS Femoral 3-in-1 Sigma™ Size 1.5	950501585	Tibial Impactor
950501533	MITKR CAS Femoral 3-in-1 Sigma™ Size 6	950501586	Femoral Introducer
950501534	MITKR CAS Femoral 4-in-1 LCS® Small	950501587	Femoral Anterior Bow Yoke
950501540	MITKR CAS Femoral 4-in-1 LCS® Large+	950501597	MI LCS® Patella Instruments - Small/Small+ Medial
950501565	MITKR CAS Patella Track Assembly	950501598	MI LCS® Patella Instruments - Small/Small+ Lateral
950501567	MITKR CAS System Screw	950501608	Tibial/Femoral RM/LL x 15 Degree Cutting Block
950501568	MI CAS NC Keel Punch Bush Size 1-1.5	950501609	MITKR CAS Patellar Wafer
950501569	MI CAS NC Keel Punch Bush Size 2-3	950501610	MITKR CAS/Specialist® 2 Tibial Adapter
950501570	MI CAS NC Keel Punch Bush Size 4-7	950501611	Quick Drill Pin - MI
950501571	MITKR CAS System Screw	950501612	MI - Patella Resection Guide - P.F.C.®
950501572	MITKR CAS Prefixation Screw	950501613	MI- Patella Resection Guide - LCS®
950501573	MITKR CAS Patella Trial - LCS® Small	129901061	Tibial Screws
950501574	MITKR CAS Patella Trial - LCS® Small+	129901030	Manual Driver
950501575	MITKR CAS Patella Trial - LCS® Standard	950501606	Tibial/Femoral LM/RL x 15 Degree Cutting Block
950501576	MITKR CAS Patella Trial - LCS® Standard+	217842005	Completion Qdrill Driver
		217842010	Completion Qdrill Pin Pack

950501599 MI LCS® Patella Instruments - Standard/Standard+ Medial
950501544 MITKR CAS 0 Degree 45 mm Cutting Block
950501546 MITKR CAS 15 Degree 45 mm Cutting Block
950501602 MI LCS® Patella Instruments - Large/Large+ Lateral
950501604 Tibial/Femoral 65 mm x 15 Degree Cutting Block
950501550 MITKR CAS Femoral 3-in-1 RP-F Size 1
950501551 MITKR CAS Femoral 3-in-1 RP-F Size 1.5
950501557 MITKR CAS Femoral 3-in-1 RP-F Size 6
950501600 MI LCS® Patella Instruments - Standard/Standard+ Lateral
950501601 MI LCS® Patella Instruments - Large/Large+ Medial
950501603 Tibial/Femoral 65 mm x 0 Degree Cutting Block
950501605 Tibial/Femoral LM/RL x 0 Degree Cutting Block
950501607 Tibial/Femoral RM/LL x 0 Degree Cutting Block
950501549 MITKR CAS Tibial/Femoral Jig Clamp
129901052 CAS Ligament Tensor Size 2 (12 Kg)
129901053 CAS Ligament Tensor Size 4 (12 Kg)
129901054 CAS Ligament Tensor Handle
129901055 CAS Tensor Spacer 5 mm
129901056 CAS Tensor Spacer 10 mm
CAS Ligament Tensor Size 2 (18 Kg)
CAS Ligament Tensor Size 4 (18 Kg)
CAS Ligament Tensor Size 2 (24 Kg)
CAS Ligament Tensor Size 4 (24 Kg)

Saw blades

950501591 MITKR CAS Saw Blade - Old Stryker
950501592 MITKR CAS Saw Blade - New Stryker
950501593 MITKR CAS Saw Blade - Old Zimmer
950501594 MITKR CAS Saw Blade - New Zimmer
950501595 MITKR CAS Saw Blade - AO/Sod
950501596 MITKR CAS Saw Blade - 3M
229910000 MITKR CAS Double Side Recip Stryker
229910001 MITKR CAS Double Side Recip Zimmer
229910002 MITKR CAS Double Side Recip Hall

INTELLIGENT  ORTHOPAEDICS™



This publication is not intended for distribution in the USA

The Ci™ logo, Ci™ and Intelligent Orthopaedics™ are trademarks of DePuy International Ltd.
Complete™, DuoFix™ and Sigma™ are trademarks and LCS® and P.F.C.® are registered trademarks of DePuy Orthopaedics, Inc.

© 2005 DePuy International Limited. All rights reserved.

Cat No: 9095-31-000 Version 1



DePuy International Ltd
St Anthony's Road
Leeds LS11 8DT
England
Tel: +44 (113) 387 7800
Fax: +44 (113) 387 7890

