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**TOTAL KNEE PROSTHESIS**  
**ANAKINE** ANATOMICAL KINEMATICS

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## CONTENTS

<b>1. Indications.....</b>	<b>4</b>
<b>2. Features.....</b>	<b>4</b>
<b>3. Surgical technique.....</b>	<b>5</b>
<b>3.1. Surgical approach.....</b>	<b>5</b>
<b>3.2. Femur preparation.....</b>	<b>5</b>
3.2.1. Preparing the intramedullary canal.....	5
3.2.2. Alignment assembly of the femur.....	6
3.2.3. Distal femoral resection.....	7
3.2.4. Femoral component measurement and positioning.....	8
3.2.5. Anterior, posterior, chamfer and trochlear femoral resection.....	10
<b>3.3. Tibia Preparation.....</b>	<b>12</b>
<b>3.3.1. Extramedular alignment.....</b>	<b>12</b>
1. Extramedular tibial guide placement and alignment.....	12
2. Tibial resection depth measuring.....	14
3. Tibial resection head fixation.....	14
<b>3.3.2. Intramedular alignment.....</b>	<b>15</b>
1. Intramedular tibial guide assembly and placement.....	15
2. Tibial resection depth measuring .....	16
3. Rotational alignment and tibial resection head fixation. ....	16
<b>3.3.3. Proximal tibial resection.....</b>	<b>17</b>
<b>3.3.4. Balancing the flexion and extension gaps.....</b>	<b>18</b>
<b>3.4. Placing and final testing of the trial components.....</b>	<b>18</b>
3.4.1. Measurement with the tibial trial tray.....	18
3.4.2. Trial femoral component placement.....	19
3.4.3. Trial reduction.....	19
3.4.4. Alignment checking.....	20
<b>3.5. Femur and tibia final preparation.....</b>	<b>21</b>
3.5.1. Tibial preparation.....	21
3.5.2. Femoral holes preparation.....	21
<b>3.6. Patella Re-Surfacing.....</b>	<b>22</b>
3.6.1. Patella preparation.....	22
3.6.2. Patella measurement.....	22
3.6.3. Patellar surface resection.....	23
3.6.4. Patellar perforation.....	23
<b>3.7. Prosthesis implantation.....</b>	<b>23</b>
3.7.1. The tibial component.....	23
3.7.2. The femoral component.....	24
3.7.3. The patella component.....	25

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*The total knee arthroplasty is indicated in those patients who have a particularly painful joint or in which it is disabled as a result of an osteo-arthritis, an arthritis caused by a trauma or a rheumatoid arthritis. The same can be said for the patients suffering constant or recurrent pain, who suffering a certain physical deterioration may need the revision of a previous arthroplasty. In the same case are different pathologies which clinical treatment is difficult and in which experience has proved that conventional techniques of reconstruction does not provide satisfactory results or in which arthrodesis is contraindicated. The knee arthroplasty should only be done when, having a progressive disability, other more conservative treatments have been used without a symptomatic improvement.*

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*The use of this total knee prosthesis is not indicated in cases in which lateral ligaments are not stabilising.*


*Although the ANAKINE total knee prosthesis has been mainly developed for primary indications, in some cases it could be implanted as a revision of a previous arthroplasty, whether it is partial or total, and for this purpose a pre-clinical evaluation of the case should be done to check if the bone, ligament and aseptic conditions allow its implantation.*

*At a tibial level, there are different modules (tibial post plug, tibial stems of different lengths and diameters) with the purpose of adapting the tibial tray to the bone conditions of the patient.*

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*A medial skin incision is made across the knee [Figure]. We get to the joint by an internal peripatellar capsulotomy. In the proximal direction, the capsulotomy goes to the medial of the quadriceps tendon. Distally, it goes down to the tibial tuberosity, in a medial direction to the patellar tendon. An eversion of the patella is done and the knee is flexed. The femoral osteophytes are removed from the intercondylar area to identify with precision the intramedullary canal. In this moment, the soft tissues are released in case of important contractions of them*

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
Figure 1.- Surgical approach

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*With the knee flexed and using the intramedullary awl (Ref.-57.35.02), open the femur entering in the intramedullary canal in the medial line of the femoral trochlea, 0,5 cm above the femoral insertion of the ACL. Cut through the cortical with the awl clockwise [Figure 2].*

*Enlarge the hole using a 9.5mm diameter drill (Ref.57.35.90), avoiding a femoral cortical wall piercing. To make the operation easy, the canal should be screwed parallel to the femoral diaphyseal axis, both in the anterior and lateral sides.*

Figura 2.-Abertura do canal medular

Figure 2.- Preparing the intramedullary canal

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*Assemble the IM reamer (Ref.-57.35.05) with the "T" handle (Ref.-57.35.04) taking away the sleeve towards the handle, **[Figure 3]** and slowly ream the intramedullary canal, cutting in a clockwise direction **[Figure 4]**. This procedure not only allows to reduce the pressure while introducing the intramedullar rod but it also facilitates its introduction. Remove the reamer and separate it from the handle.*

*Assemble the intramedullar rod (Ref.-57.35.06) onto handle and slowly insert it with a torsion movement to eliminate the canal content. Release the "T" handle leaving the rod in the canal **[Figure 5]**.*

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*Determine the appropriate angle in valgus of the patient femoral diaphysis in a pre-surgical way with the appropriate template, over an extra-long xray of the leg in extension and select mentioned angle in the femoral alignment angle guide (Ref.-57.35.08) pressing the screw until it is fixed **[Figure 6]**. The body of the angle guide can be adjusted with increases of 1 degree, from 11 Left to 11 Right, which allows the correct placing of the distal femoral resection guide (Ref.-57.35.09). Special care should be taken when keeping the blocking screw tight in the correct position*

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*Place the alignment guide in the IM bar and slide it until touching the distal femoral condyles. In cases of severe condylar erosion, only the least eroded condyle should touch the face of the guide. You should not attempt to alter the guide from its pre-set angle, because it can result in a misaligned femoral component [Figure 7].*

*Rotate the distal femoral alignment guide so the posterior margins are parallel to the femoral posterior condyles. Once it is correctly in its position, place the 3.2 mm diameter anchoring pins to fix it (Ref.-57.35.88/89). If the rotational alignment is not correct, this can negatively affect the varus-valgus alignment.*

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*Insert the graduated support with the pivots through the holes in the femoral alignment angle guide (Ref.-57.35.08) [Figure 8]. Assemble the external alignment tower (Ref.-57.35.16) to check external alignment [Figure 9]. The alignment rod (Ref.57.35.14) should face the centre of the femoral head. Once the external alignment is checked, remove the alignment tower and the rod*

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*Slide the distal femoral resection guide (Ref.-57.35.09) over the support in the angle guide by its connexion system.*

The resection guide, is oriented so the marks are on the viewed surface [Figure 10].

This support has divisions from 3mm to 25mm allowing the guide to be placed at the required resection level; at this point, the screw should be tight to secure the position in which the femoral resection guide is fixed. To resect the same thickness of the condyles in the femoral component, place the mark 0 at the 9mm level (value rounded with a circle) [Figure 11].

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It is important to re-establish the position of the trochlear sulcus, which can be achieved by touching with the reference finger (ref.-57.35.10) inside the intercondylar groove [Figure 12].

Hold the support and the resection guide on the anterior cortical. Introduce the anchoring pins in the holes, marked with 0 and make them go forward to the anterior cortical [Figure 12].

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After introducing the pins, it is still possible to correct the resection level, replacing the guide in the distal or proximal direction, with the pins, in increments of 2mm or using the cutting slots marked with +4 or +8. Table 1 shows us the different possible corrections. For a resection equal to the selected level in the graduated support, you should use the holes and slot of level 0 (value inside the square).

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*The cutting block has some oblique lateral holes to, once the anchoring pins are placed, avoid the migration of the guide when making the cut. It also has proximal holes for increased stability. Once the resection guide is safe in its position, you can disassemble, following this order, the IM rod, the angle guide and the support. To remove the pins used during surgery, use the nail extractor (Ref.-57.35.13). To facilitate the insertion of the anchoring pins you can use the cancellous bone drill (Ref.-57.20.20).*

*Resect the distal femoral bone using an oscillating saw blade through the cutting guide slots **[Figure 13]**. After making the cut, you should work the saw from back to front in order to finish a flat surface at the level of the slot.*

⚠ *You must use 1,2mm thick saw blades in order to achieve a precise cut with the guides. If the surgeon decides to cut over the block surface, 4mm less of bone will be cut, and that should do the corresponding adjustment.*

*You should level any irregularity in the osteotomy surface with the bone file (Ref.-57.20.94).*

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The antero/posterior, medial/lateral and rotational position is determined by the placement of the femoral drill/sizing guide (Ref.-57.35.18) [Figure 14].

Place the guide at the level of the distal femoral cut surface. The guide feet should sit firm touching the posterior condyles. If there is a posterior condyle defect, the sizing guide should be turned to a position compensating mentioned defect. Position the stylus point over the anterior cortex, and indicate the level at which the anterior femoral resection will emerge. Block the position determined with the screw, and check that the stylus is still placed in the anterior cortex, without protrusion. The estimated size of the femoral component is indicated on the femoral drill guide [Figure 14, 15]. If the reading falls between two sizes select the larger option. The x-ray templates (Ref.-57.35.07) are provided to determine optimal size selection. Once the appropriate size is chosen, you should centre the guide in a medial-lateral position. As an option, we can use the lateral finger, sliding it up to the selected size with the sizing guide, and holding it against the lateral cortical of the femur, making use of the 3.5 mm screwdriver (Ref.-66.00.02). In this way, the femoral component is implanted in the lateral side, at the level of the lateral cortical [Figure 14].




Figura 14.- Posicionamento do guia de perfuração / medidor femoral  
Figure 14.- Femoral drill/ sizing guide placement

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To establish a correct rotation of the femoral component it helps to use the clinical epicondylar axis. You should select a module of femoral drill for a neutral, 3 or 4 degrees external rotation (Ref.-57.35.19/21). Insert the module in the appropriate side (right or left) and make the holes using the femoral drill with stop dia Ø4mm (Ref.57.35.17) [Figure 16], after the first hole has been prepared, insert the Tommy bar (Ref.-57.35.15) to keep proper alignment for the next hole. The instrumental provides a reference finger (Ref.-57.35.10) which shows the end of the femoral anterior cut.

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Figure 16.- Drilling through a femoral drilling module

### 3.2.5. Ressecção femoral anterior, posterior, chanfradura e tróclear / *Anterior, posterior, chamfer and trochlear femoral resection*

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*Once the femoral drill/sizing guide is taken away and the correct size of the component is determined, you should select the appropriate size resection guide (Ref.-57.35.22/26) and introduce its pins into the prepared holes in the distal femur [Figure 17]. Introduce two anchoring pins in the lateral holes to provide more stability. Care should be taken to minimize translation forces on the cutting guide which tend to enlarge the alignment holes. Make the anterior, posterior and chamfer cuts with an oscillating saw. Place the saw blade through the slot prior to activating it. First resect the anterior and posterior condyles, then complete the chamfer cuts.*

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Surgical technique ANAKINE TOTAL KNEE PROSTHESIS  
 Use the resection guide handles (Ref.-57.35.27) to make it stable when making the cuts. The handles have a superior bending to avoid grabbing the soft tissues during the insertion, but if necessary you can remove them to make the resection easier.

Remove the cutting guide and check the cut surfaces. If it is necessary, level any surface irregularity with the bone file. Then, secure the intercondylar chisel Trochlear guide (Ref.-57.35.42/46) with the size previously determined, to the femur distal surface inserting the guide pins in the holes made [Figure 18]. The guide should be well located on the anterior cortical and the distal surface. For this purpose you should use both resection guide handles (Ref.-57.35.27), making a pressure on these

surfaces. Introduce the anchoring pins in to the oblique holes of the guide, in the anterior face. Introduce the intercondylar chisel of the appropriate size (Ref.57.35.52/53) in to the opening marked as "Trochlear" in the intercondylar resection guide. The chisel is oriented so the marks appear on its surface and it advances until its stop..

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*Anakine total knee system offers the option of referencing the tibial mechanical axis with an IM tibial alignment guide or an external tibial alignment guide. In patients whose tibial medullary canal is severely bowed, the IM alignment approach may tend towards malalignment. In these cases, use of the external tibial alignment guide may be indicated. Preoperative x-ray evaluation may assist in this selection.*

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*Place and stabilise the knee in complete flexion applying distraction to the tibia in anterior direction. Choose the appropriate tibial posterior slope module (Ref.-57.36.01/02). There are two types of interchangeable modules provided, one identified as 0 DEG that provides a neutral resection and another identified as 3 DEG which*

provides a posterior slope of 3° in the tibia sagittal plane. Place the selected module in the extramedullar guide (Ref.-57.36.04) and over it the tibial resection head (Ref.-57.35.99/57.36.00) fixing the position with the screw [Figure 22]. The tibial resection head is anatomical and it is available in two versions, depending on the operated knee being the left or right (Ref.-57.35.99 / 57.36.00).

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Figure 22.- Lateral plane of the extramedullar tibial resection / Traverse plane; anatomical tibial resection head

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Figura 23.- Plano frontal do ressecor tibial extramedular

Figure 23.- Front show of the extramedullar tibial resector

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Figure 24.- Medial-lateral alignment of the ankle joint clamp

Place the EM tibial resector in the tibial front part and fix it surrounding the ankle with the ankle clamp, placing it in the immediate proximal area to the malleolus. Level the platform at the proximal plateau. Move the lower part of the assembly in the anterior-posterior direction to align it parallel to the tibial axis [Figures 22 and 23].

The medial-lateral alignment is nearly parallel to the tibial axis, place the lower assembly slightly into medial

position [Figure 24]. In the cases where the resection head has moved, the lower assembly should be adjusted, moving it in the same direction. Once the extramedular tibial guide is placed, it can be stabilized with anchoring pins.

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*Assemble the Stylus (Ref.-57.36.03) in the tibial resection head by the holes done for such purpose, so the tip marked as “most normal” or that marked as “most defective” go over the tibial plateau [Figure 25].*

*If the surgeon decides to measure the resection depth from the tibial plateau most integral or normal, the “most normal” side of the stylus should be placed in the centre of the less damaged condyle (less depressed). This level provides 10 mm of bone resection below the point marked by the stylus.*

*If the resection depth measurement is done from the most defective (most depressed) side, the “most defective” tip of the stylus should be placed in the centre of the most defective tibial side. The bone resection level is 2mm below the point marked by the stylus.*

*The extramedular alignment guide has a screw that allows a millimetrical proximal-distal adjustment of the resection head [Figure 26].*

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
*Once the correct placing of the extramedular alignment guide is guaranteed, introduce the anchoring pins into the tibia through the tibial resection head. Introduce them into the holes marked with . This way, the resection will be of 10mm from the "most normal" or 2mm from the "most defective" regarding the measurement made by the Stylus. The resection head could be moved 2mm proximally or distally, in order to increase or decrease the tibial resection. The block has some oblique lateral holes that avoid the guide migration when making the cut [Figure 27]. Once the head is fixed, remove the extramedular alignment guide unlocking the anterior screw.*

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*Place and make the knee stable in flexion applying distraction to the tibia in anterior direction. To open the IM tibial canal you can follow the same steps used in the femur.*

*Locate the entry point with the maximum precision, finding the proximal extension of the IM canal centre in the pre-surgery x-ray. This point is normally placed between the tibial spines in the medial-lateral side, just anterior to them in the anterior-posterior direction.*

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Figure 28.- Anatomic tibial resection head

*With the IM in the canal, slide the tibial resector IM (Ref.-57.36.50) through the bar and place it at the level of the tibial surface. For the guide to be placed easy in the tibia proximal surface, you can resect the tibial spine at the highest level of the less defective tibial plate.*

*Similar to the extramedular alignment, there are two posterior slope modules available (Ref.-57.36.51/52), one of 0º and another of 3º respectively, which are placed in the tibial alignment guide and are used as a base for the tibial resection head. Hold the resection head on the anterior cortical fixing the slope module in the guide. The tibial*

*resection head is anatomical and available in two versions, depending on the operated knee being the left or the right (Ref.-57.35.99 / 57.36.00) [Figure 28].*

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*Assemble the Stylus (Ref.-57.36.03) in the tibial resection head through the holes done for such purpose, so the tip marked as “most normal” or that marked as “most defective” goes over the tibial plateau.*

*If the surgeon decides to measure the resection depth from THE RESECTION the tibial plateau most integral or normal, the “most normal” side of the stylus should be placed in the centre of the less damaged condyle (less eroded). This level provides 10 mm of bone resection below the point marked by the stylus.*

*If the resection depth measurement is done from the most defective (most eroded) side, the “most defective” tip of the stylus should be placed in the centre of the most defective tibial side. The bone resection level is 2mm below the point marked by the stylus [Figure 29].*

*Adjust the position of the resection head in an anterior-posterior and proximal-distal direction and secure the screws placed in the IM guide. Block the screw of the resection head [Figure 30].*

Figura 29.- Colocação do ressector tibial intramedular  
*Figure 29.- Placing of the IM tibial resector*

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*Determine the medial-lateral alignment connecting the external alignment tower (Ref.-57.35.16). Pass the extramedular alignment rod (Ref.-57.35.14) through the most appropriate hole of the tower. The rod should be aligned with the ankle joint centre.*

*Once the resection head placement is correct, insert the anchoring pins in the tibia through the tibial resection head. Introduce them in to the holes marked with a . This way, the resection will be of 10mm from the “most normal” or 2mm from the “most defective” regarding the measurement made by the Stylus. The resection head could be moved 2mm proximally or 2mm distally, in order to increase or decrease the tibial resection. The block has some oblique lateral holes that avoid the guide migration when making the cut. Once it is fixed, remove the intramedular alignment guide unlocking the anterior screw of the head and extract the posterior slope module.*

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*Cut the proximal tibia using an oscillating saw inserted through the cutting slot of the resection head **[Figure 31]**. After making the cut you should work the saw from back to front in order to finish the surface at the level of the slot.  You must use a 1.2mm thick saw to achieve a precise cut with the guides.*

*The surgeon should pay special attention during the tibial resection to avoid damaging the collateral and posterior, nervous and vascular structures. The tibial fixation of the PCL must be kept full or partially integrated when using a PCL preservation prosthesis (CR).*

*Level any irregularity of the resected surface with the saw or with a bone file (Ref.-57.20.94).*

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#### **BALANCING THE FLEXION AND EXTENSION GAPS**

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*Assemble the appropriate flexion-extension gauge (Ref.-57.36.32/35) to the tibial handle (Ref.-57.36.46) and check the ligament tension in flexion and extension [Figure 32].*

*The main principles in the balancing of the flexion and extension gap are:*

- 1 *The amount of distal femur resected determines the soft tissue tension in extension.*
- 2 *The amount of posterior femoral condylar bone determines the soft tissue tension in flexion.*

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Figure 32- Balance of the extension gap

3. *The amount of tibial bone resected determines the soft tissue tension in both flexion and extension. Correct alignment should be achieved by the bone resections. However, if it is still necessary to further leg alignment by either a medial or lateral collateral ligament release, care should be taken and the flexion/extension gaps checked again.*

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*With the tibia in anterior subluxation position, select the tibial trial tray (Ref.-57.36.36/40) that offers better coverage of the prepared surface, assemble the tibial handle (Ref.-57.36.46) and place it over the resected tibial surface. The trial tray should provide an optimal coverage of the tibia in the medial-lateral and anterior-posterior edges [Figure 33]. Table 3 shows the possibilities of interchanging between sizes of tibial trays and inserts.*

Bandeja TRAY	Meseta INSERT	1	2	3	4	5
1		si YES	si YES	no	no	no
2		si YES	si YES	si YES	no	no
3		si YES	si YES	si YES	si YES	no
4		si YES	si YES	si YES	si YES	si YES
5		si YES	si YES	si YES	si YES	si YES

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*Insert the femoral impactor-extractor (Ref.-57.35.41) to the universal handle (Ref.-57.35.39) by its fast connexion mechanism. Assemble the femoral trial to the impactor-extractor and place it over the femoral surface. The anterior edges are taken forward at the same time, parallel to the distal femoral cut. To finish the impact, use the femoral impactor (Ref.-57.35.38) assembled to the universal handle.*

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*Introduce the tibial trial insert (CR or PS depending on needs Ref.-57.36.12/31 and Ref.57.36.41/45),  the size must coincide with that of the femoral trial previously selected. Once the trials are place in their position extent the knee slowly evaluating the medial and lateral stability and the alignment in the anterior-posterior and medial-lateral planes. If the surgeon finds that there is instability, the tibial trial insert should be changed with a thicker one, making a new reduction. Choose the tibial thickness that provides better stability in the flexion-extension, allowing a complete extension.*

*Adjust the rotational alignment of the tibial trial tray with the knee in extension. Make use of the tibial handle to rotate the tibial trial assembly tray-insert and make them coincident with the femoral trial component. Mark with a electrocautery device the appropriate position over the anterior tibial cortical [Figure 34].*

*With the knee in complete flexion and with the tibia in anterior subluxation position, remove the tibial trial insert. Test with the trial tray assembled in the tibial handle the appropriate rotational alignment taking as reference the electrocautery marks. Fix the tray using the headed tibial anchoring pins (Ref.-57.36.53/54) and connect the trial insert previously selected on the tibial trial tray.*

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*With the knee in complete extension, assemble the two alignment rods (Ref.-57.35.14) and introduce them through the appropriate hole. Assess the alignment [Figure 35]. The correct rotation of the tibial component should be determined in agreement with its position with the femoral component. When the alignment is correct, the rod passes through the mechanical axis of the hip, knee and ankle. Once the alignment is checked, flex the knee completely and bring the tibia to an anterior subluxation position again. Remove the assembly leaving in position the trial tray, the handle and the pins.*

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*Select the appropriate tibial punch-drill guide (Ref.-57.36.05/06) and place it in the space prepared for it on the tibial trial tray. Fix the screw of the guide **[Figure 36]** with the 3.5 mm screwdriver (Ref.-66.00.02). Pass the tibial drill (Ref.-57.36.09/10) of the appropriate size adjusting the depth stop in the mark “PLUG, “STEM 25” or “STEM 50” **[Figure 37]**. Assemble the tibial punch (Ref.-57.36.07/10) to the universal handle (Ref.-57.35.39) by its fast connexion mechanism and introduce it through the guide pushing it to the stop. Disassemble the handle and assemble the slap hammer (Ref.-57.35.40) to remove the punch.*

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*Check the medial-lateral position of the femoral trial component and prepare the holes for the femoral pegs, advancing the femoral drill with stop dia 6,3mm (Ref.-57.35.11) or the femoral punch (Ref.-57.35.12) through the holes of the femoral trial **[Figure 38]**.*

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*Remove in a circle the osteophytes in the patella peripheral and the synovial insertion, respecting the insertion area of the quadriceps and the patella ligament. Release enough soft tissue in the pre-patella bursa in order to place the patellar caliper in the anterior cortical.*

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*Measure the original thickness of the patella by placing the patella caliper tips (Ref.57.37.01) over the sides of the joint surface and anterior cortex [Figure 39]. The biggest sagittal thickness is on the medial edge. The average range is between 20 and 30 mm. The average value is taken from the line marked with “OUT”.*

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*Select the patella template (Ref.-57.37.02/03) that covers the joint surface in the most appropriate way, so it does not overlap out. Place the handle in the medial side of the everted patella. You should not use a small patella (32x27) with the longest femoral component (size 5).*

*Take the patella original thickness measured with the patella calliper and adjust the level of resection with the gauge placed in the resection guide (Ref.-57.37.08) [Figure 40]. In case of using the sizes of 38x33 or 40x36, subtract 2mm of mentioned original thickness.*

*Hold the guide parallel bars over the anterior cortex and place the edge dent side on the upper and lower margins of the articular joint surface. Close the handles to hold the patella firmly.*

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*Make the resection with the oscillating saw, holding it at the level of the cutting surface. The resection should be made through the slots in the guide head [Figure 41]. Then, remove the guide and check the residual thickness with the patella caliper.*

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*Place the appropriate patela template (Ref.-57.37.02/03) over the cut surface placing the handle in the medial side of the everted patella. Press the template firmly on the resected surface and make the holes with the patella drill with stop (Ref.-57.37.00) [Figure 42]. To check the correct position and implant bone coverage, you can use the patella trials (Ref.-57.37.04/07) in this moment or during the reduction trial.*

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*The tibial tray has a pre-assembled closing plug that could be taken away in case of needing tibial stems. Make*

use of the tibial stem wrench (Ref.-57.36.47) to release the plug and fix the stems. Clean the resected surface and apply some cement covering the underside of the tibial tray with all its surfaces. Assemble the universal handle to the tibial impactorextractor (Ref.-57.36.48) and the assembled universal handle to the definitive tibial tray. Carefully insert the tray avoiding a bad rotation. As an alternative to the tibial impactorextractor, you can use the tibial impactor (Ref.-57.36.49). Once the tray is fully seated, slap the upper part of the universal handle until fixing it **[Figure 43]**. Remove the excessive cement of the edges, beginning from the back and moving towards the sides and the front part.

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Slide the definitive tibial insert (CR or PS, depending) over the tray in the anterior-posterior direction to the stop **[Figure 44]**. **⚠ The size of the insert must be the same as the femoral component.** Fix the two components with the locking screw using the 3.5 mm screwdriver (Ref.-66.00.02) **[Figure 45]**.

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*Clean the resected surface.*

*If a cemented femoral component is used, apply cement both on the bone surface and in the femoral component. Assemble the femoral component (CR or PS as required) to the femoral impactor-extractor (Ref.57.35.41). The holding edges are taken forward at the same time, parallel to the distal surfaces and protecting the surfaces carefully set up until completely seating the interior surfaces. Release the device and complete the proceeding with the femoral impactor (Ref.-57.35.38) and the slap hammer (Ref.-57.35.40) **[Figure 46]**. Clean the rests of cement if necessary.*

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*Clean carefully the patella resected surface. The patellar implantation can be cemented at the same time as any of the other components. Apply cement on the cut surface. Apply the curved surface of the patellar clamp (Ref.-57.37.09) over the implant, making pressure **[Figure 47]**. Take away the excessive cement from the implant edges, and hold the compression with the clamp, until the cement is dry.*

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