




# CORAIL<sup>®</sup>

---

HIP SYSTEM

---

Surgical Technique



## *Intelligent Surgery*

*Intelligent hip surgery is an approach to patient treatment that places equal importance on:*

*Optimising function*

*Maximising survivorship*

*Accelerating recovery*

*Its success is founded on leadership in the development of:*

*High performance bearings*

*Clinically proven implants*

*Responsible and effective MI techniques*

“In use worldwide; FDA approved in 1996; recognised Level 10A in 2004 by the Orthopaedic Data Evaluation Panel in the UK: the Corail® stem has now become a gold standard among primary stems. The Corail® stem’s philosophy is based on simple principles: primary mechanical stability, secondary biological integration, bone preservation and harmonious stress transfer. The design, unchanged since 1986, gives the primary mechanical stability. The hydroxyapatite coating allows secondary biological integration. The combination of the design and the HA coating of the Corail® stem has proven to work perfectly.<sup>1-3</sup> The surgical technique is simple and allows bone preservation as we are looking for an “optimum filling” and not a close cortical contact with the implant. The restoration of bone stock occurs with the creation of newly formed bone all around the stem thanks to the effect of both the design and the hydroxyapatite. The compaction broaching surgical technique is reproducible and straightforward. We do not see any long-term radiographic changes.”

ARTRO Group  
Corail® Design Surgeon Team  
Clinique d’Argonay  
International Visitation Centre Corail®  
Annecy, France




Pre-op. A / P view.



Post-op. Corail® stem with 36 mm Ceramax™ Head and Pinnacle™ Cup System.

After twelve years of constant pain, Johan was unable to sleep through the night and facing the fact that he could no longer run his business as before, Johan's quality of life at 44 years old was extremely poor. Life was no longer fun. Now, just seven months after his Corail® Hip surgery, he is back enjoying life with friends, skiing and ice climbing.

# CLINICALLY PROVEN HIP



95.1%

Survivorship in 5,130 cases\* at 15 years. The Norwegian Arthroplasty Register 1987-2004, 2005 <sup>1</sup>

98.3%

Survivorship in 2,956 cases at 10 and 14 years. Vidalain JP. Artro Group., 1998 <sup>2</sup>

98.9%

Survivorship in 100 consecutive cases at 8 years. Røkkum M., J. Bone and Joint Surg., 1999 <sup>3</sup>

\* Cohort size on Norwegian Hip Arthroplasty Register.

# PRE-OPERATIVE PLANNING


## Pre-operative planning

The Corail® Hip System provides pre-operative templates at three different magnifications (100%, 115% and 120%). These are placed over the A/P and M/L radiographs to help determine the implant size in order to restore the patient's natural anatomy. The objective is not for the implant to fill the femur in close cortical contact, but to sit the stem in the compacted cancellous bone. There should be at least a 1 mm distance between the stem and the cortical bone. The pre-operative templating will also indicate the level of neck resection.

## Surgical approach


The Corail® stem can be used with any surgical approach using a conventional or a reduced incision.

# FEMORAL NECK RESECTION



The angle of resection should be 45°. The neck resection guide should be used to determine the level of the femoral neck resection in conjunction with pre-operative templating. If the resection is too high, it may result in a varus positioned stem. *Note: the osteotomy can be performed in one or two planes depending on the surgeon's preference.*


# PROXIMAL CANCELLOUS BONE COMPACTION



Enter the femoral canal as laterally as possible with a chisel to avoid varus positioning. Use the bone tamp to compact the cancellous bone proximally. This is an important step as the philosophy of the Corail® stem is based on bone preservation. Please refer to the Pinnacle™ Surgical Technique for full details with regards to the acetabulum preparation (cat no: 9068-80-050).




# FEMORAL CANAL PREPARATION




Begin with the smallest broach attached to the broach handle and increase the size of broach one at a time until axial and rotational stability is achieved. Stop broaching when it is felt that the broach is stable longitudinally and rotationally. Careful pre-operative planning is key to help selection of the final broach size. Anteversion is automatically set during the broaching with the flared broaches.

# CALCAR REAMING




Leave the last broach in place and use the calcar mill to achieve a flat resection surface. If a collared stem is used, the calcar reaming should allow an optimised fit of the collar on the calcar.

# TRIAL REDUCTION



With the last broach in situ, attach the appropriate trial neck and trial head. Reduce the hip and assess what adjustments, if any, are required to ensure stability through a full range of motion. Remove the femoral head, neck trial and broach. Do not irrigate or dry the femoral canal. This will help to preserve the compacted cancellous bone quality and encourage osteointegration of the stem.

# FEMORAL COMPONENT INSERTION



The impactor is used to impact the definitive stem for the last few centimetres.

A 3D anatomical illustration showing the femoral component insertion process. A femoral head and neck are shown in a light tan color. A pink femoral head is seated within a metal acetabular cup. A metal femoral neck is being inserted into the femoral canal. A metal impactor is positioned over the neck, with its handle extending upwards. A line points from the text to the impactor. Another line points from the text to the junction of the HA coating and the milled femoral neck.


The HA coating should sit level with the milled femoral neck.

Engage the Corail® stem (that has the same size as the last used broach) in the femoral canal by hand and finish its introduction with the impactor for the last few centimetres.

*Note: The stem is 0.31 mm thicker than the broach to allow the necessary press-fit.*


# ADDITION OF BONE GRAFT

Cancellous bone graft is packed around the proximal part of the prosthesis to encourage bone on-growth.



Once the Corail® stem is seated, cancellous bone from the resected femoral head is added around the proximal part of the stem using the bone tamp to seal the femoral canal and to reduce the time for osteointegration which provides definitive stability.

# FEMORAL HEAD IMPACTION



Clean the stem taper carefully to remove any particulate debris. Place the femoral head onto the taper and lightly tap it (especially if a ceramic head is used) using the head impactor. Ensure bearing surfaces are clean, and finally reduce the hip.


# ORDERING INFORMATION

## Corail® Implants

All the Corail® stems have the Articul/eze™ neck, characterised by a thin antero-posterior dimension, a polished surface and a 12/14 mini-taper.

### Corail® Standard Stem (Collarless)

3L92507	Corail® Size 8
3L92509	Corail® Size 9
3L92510	Corail® Size 10
3L92511	Corail® Size 11
3L92512	Corail® Size 12
3L92513	Corail® Size 13
3L92514	Corail® Size 14
3L92515	Corail® Size 15
3L92516	Corail® Size 16
3L92518	Corail® Size 18
3L92520	Corail® Size 20



### Corail® Standard Stem (Collared)

3L92498	Corail® Size 8
3L92499	Corail® Size 9
3L92500	Corail® Size 10
3L92501	Corail® Size 11
3L92502	Corail® Size 12
3L92503	Corail® Size 13
3L92504	Corail® Size 14
3L92505	Corail® Size 15
3L92506	Corail® Size 16
3L92508	Corail® Size 18
3L92521	Corail® Size 20



### Corail® Coxa Vara Lateralised Stem (Collared)

3L93709	Corail® Size 9
3L93710	Corail® Size 10
3L93711	Corail® Size 11
3L93712	Corail® Size 12
3L93713	Corail® Size 13
3L93714	Corail® Size 14
3L93715	Corail® Size 15
3L93716	Corail® Size 16
3L93718	Corail® Size 18
3L93720	Corail® Size 20



### Corail® High Offset Lateralised Stem (Collarless)

L20309	Corail® Size 9
L20310	Corail® Size 10
L20311	Corail® Size 11
L20312	Corail® Size 12
L20313	Corail® Size 13
L20314	Corail® Size 14
L20315	Corail® Size 15
L20316	Corail® Size 16
L20318	Corail® Size 18
L20320	Corail® Size 20



## STANDARD - COLLARLESS / COLLARED

Size	Stem Length (mm) (A)	Stem Length (mm) (B)	Offset (mm) (C)	Neck Length (mm) (D)	Neck Shaft Angle (E)
8	115	95	38.0	38.5	135°
9	130	110	38.5	38.5	135°
10	140	120	39.5	38.5	135°
11	145	125	40.0	38.5	135°
12	150	130	41.0	38.5	135°
13	155	135	41.5	38.5	135°
14	160	140	42.0	38.5	135°
15	165	145	43.0	38.5	135°
16	170	150	43.5	38.5	135°
18	180	160	44.5	38.5	135°
20	190	170	45.5	38.5	135°

## HIGH OFFSET - COLLARLESS

Size	Stem Length (mm) (A)	Stem Length (mm) (B)	Offset (mm) (C)	Neck Length (mm) (D)	Neck Shaft Angle (E)
9	130	110	45.5	43.2	135°
10	140	120	46.5	43.2	135°
11	145	125	47.0	43.2	135°
12	150	130	48.0	43.2	135°
13	155	135	48.5	43.2	135°
14	160	140	49.0	43.2	135°
15	165	145	50.0	43.2	135°
16	170	150	50.5	43.2	135°
18	180	160	51.5	43.2	135°
20	190	170	52.5	43.2	135°

## COXA VARA - COLLARED

Size	Stem Length (mm) (A)	Stem Length (mm) (B)	Offset (mm) (C)	Neck Length (mm) (D)	Neck Shaft Angle (E)
9	130	110	45.5	40,3	125°
10	140	120	46.5	40,3	125°
11	145	125	47.0	40,3	125°
12	150	130	48.0	40,3	125°
13	155	135	48.5	40,3	125°
14	160	140	49.0	40,3	125°
15	165	145	50.0	40,3	125°
16	170	150	50.5	40,3	125°
18	180	160	51.5	40,3	125°
20	190	170	52.5	40,3	125°


Note: All measurements are based on a 28 mm +5.0 Articul/eze® head, which is the middle length of non-skirted femoral heads

Corail® Instrumentation

CONTK2 Sterilisation Case



L20500 Base Aluminium Basket



L20504 Top Basket



L20503 Superior Thermoformed Tray



9653-68-000 Anteversion Axis



L20440 Neck Resection Guide



L93606 Bone Tamp



L93205 Bone Impactor



2002-31-000 Osteotome



9522-11-500 Curved Broach Handle



9522-10-500F Straight Broach Handle



9522-12-500F Extra Curved Broach Handle





L20502 Middle Thermoformed Tray



- L20408 Broach 8
- L20409 Broach 9
- L20410 Broach 10
- L20411 Broach 11
- L20412 Broach 12
- L20413 Broach 13
- L20414 Broach 14
- L20415 Broach 15
- L20416 Broach 16
- L20418 Broach 18
- L20420 Broach 20



L20431 Standard Neck Segment



L20432 Lateralised Neck Segment (Coxa Vara)



L20433 High Offset Neck Segment



L20501 Inferior Thermoformed Tray



- 2570-04-200 Calcar Mill Large
- 2570-04-100 Calcar Mill Small



- 2570-05-000 Positioner



- 2570-05-100 Stem Impactor



- 2001-65-000 Head Impactor



#### Pre-operative Templates

##### Pre-operative Templates (100%)

- CALQ861 Corail® Standard (28 mm, 32 mm, 36 mm, ASR™ XL Heads)
- CALQ862 Corail® Coxa Vara (28 mm, 32 mm, 36 mm, ASR™ XL Heads)
- CALQ863 Corail® High Offset (28 mm, 32 mm, 36 mm, ASR™ XL Heads)

##### Pre-operative Templates (120%)

- CALQ864 Corail® Standard (28 mm, 32 mm, 36 mm, ASR™ XL Heads)
- CALQ865 Corail® Coxa Vara (28 mm, 32 mm, 36 mm, ASR™ XL Heads)
- CALQ866 Corail® High Offset (28 mm, 32 mm, 36 mm, ASR™ XL Heads)

##### Pre-operative Templates (115%)

- CALQ858 Corail® Standard (28 mm, 32 mm, 36 mm, ASR™ XL Heads)
- CALQ859 Corail® Coxa Vara (28 mm, 32 mm, 36 mm, ASR™ XL Heads)
- CALQ860 Corail® High Offset (28 mm, 32 mm, 36 mm, ASR™ XL Heads)

#### Digital Templates

The availability of digital templates depends on DePuy International's agreement with the vendors. Please contact DePuy International for more information.





## References:

1. The Norwegian Arthroplasty Register 1987-2004, Prospective Studies of Hip and Knee Prostheses. <http://www.haukeland.no/nrl/>, AAOS, 2005.
2. Vidalain JP Corail® Stem Long-Term Results Based upon the 15-Years ARTRO Group Experience. Fifteen Years of Clinical Experience with Hydroxyapatite Coatings in Joint Arthroplasty, Ed. Springer, 217-224, 2004.
3. Røkkum M, Brandt M, Bye K, Hetland KR, Waage S, Reigstad A. Polyethylene Wear, Osteolysis and Acetabular Loosening with an HA Coated Hip Prosthesis. J. Bone and Joint Surg. 81-B, No4, 1999.

This publication is not intended for distribution in the USA

Corail® is a registered trademark of DePuy (Ireland) Ltd.

ASR™ is a trademark of DePuy International Ltd.

Ceramax™ and Pinnacle™ are trademarks and Articul/eze® is a registered trademark of DePuy Orthopaedics, Inc.

© 2006 DePuy International Limited. All rights reserved.

Cat No: 9066-35-001 version 1



**DePuy International Ltd**  
St Anthony's Road  
Leeds LS11 8DT  
England  
Tel: +44 (113) 387 7800  
Fax: +44 (113) 387 7890

**DePuy France S.A.S.**  
7 Allée Irène Joliot Curie  
69800 Saint Priest  
France  
Tél : +33 (0) 472 79 27 27  
Fax : +33 (0) 472 79 28 28

