

**For the Educational Use of Howmedica Osteonics
Employees Only. Not for Reproduction or Distribution.**



Surgical Technique

**The
Minimally
Invasive
Uni Knee
System**



stryker
Howmedica
OSTEONICS



Stryker® Howmedica Osteonics would like to thank the following surgeons for their help in developing and evaluating the EIUS™ Unicondylar Knee System and this surgical protocol:

Anthony DeBari, M.D.

Robert Gorab, M.D.

John Kilgore, M.D.

Richard Levitt, M.D.

Paul Lotke, M.D.

Craig McAllister, M.D.

Thomas Mallory, M.D.

Terence Matthews, M.D.

Daniel Murphy, M.D.

Matthew Phillips, M.D.

Gary Poehling, M.D.

Michael Reilly, M.D.

Todd Schmidt, M.D.

Keith Schroeder, M.D.

Enno Steinheisser

CONTENTS	PAGE NO
Introduction	
System Features/Benefits	2
Surgical Technique	
Clinical Indications	4
Basic Bone Cuts	4
Preparation	5
Surgical Steps	
Tibial Alignment	6
Tibial Resection	7
Femoral Preparation	11
Trial Reduction	14
Tibial Keel Preparation	15
Final Implantation	16
Key to Instruments/Implants	
Tray Layout	17
Femoral/Tibial Components	19

Bone Sparing Implant Design

EIUS[™] Femoral Component

- Bone sparing, resurfacing design
- Anatomical bend and contoured anterior lip designed to promote optimal condylar fit and help prevent patella impingement

EIUS[™] Tibial Component

- Designed to accommodate +/- 10° varus/valgus malalignment
- 1 up/1 down size interchangeability with femoral component



Precision Instrumentation

Tibial Preparation

- "Tibia-first" approach designed to promote accurate femoro-tibial alignment
- Low-profile, extramedullary guides designed to facilitate a minimally invasive procedure

Femoral Preparation



- Unique, patent-pending distal burring template and specialized bur designed for accurate and reproducible femoral resurfacing
- Extramedullary guides designed for true, minimally invasive technique with accurate alignment

Minimally Invasive with Confidence!

The EIUS™ Uni Knee combines a bone sparing implant and precision instrumentation to create a minimally invasive system with the ease of use of traditional "open" procedures.

Advantages include:

- Potential for shorter hospital stays and quicker rehabilitation^{1,2}
- A small, 8-10 cm incision
- No eversion of the patella
- No reaming of the intramedullary canal



¹ Romagnoli, Sergio. The Adventures of Unicompartmental Knee Replacement. Proceedings from Current Concepts in Joint Replacement. Spring 2001, Paper 33.

² Murry, David. Unicompartmental Knee Replacement: Now or Never? Orthopedics. September 2000, Vol. 23, page 980.

Appropriate indications are critical for a successful outcome in unicompartmental knee arthroplasty.

Recommended indications include:

- Osteoarthritis (medial or lateral)
- Unicompartmental post traumatic osteoarthritis
- Avascular Necrosis - medial femoral or tibial condyle
- Ahlback I-III

Possible contraindications include:

- Tricompartmental disease
- Multidirectional/ACL instability
- Severe tibial bone loss/deformity (over 15 degrees varus)
- Significant flexion contracture (over 12 degrees)
- Inflammatory arthritis
- Symptomatic patello-femoral disease or tracking abnormalities

The use of weight-bearing radiographs to determine the arthritic condition of the knee is recommended. These include the use of A/P, lateral and skyline views. Long-leg standing films can also be used for determination of the mechanical axis. Diagnostic arthroscopy of the opposite compartment and patello-femoral joint is optional as a pre-operative assessment.

Six Basic Bone Cuts

1 Sagittal tibial



2 Transverse tibial



3 Posterior femoral



4 Femoral peg and fin




5 Distal femoral burring




6 Tibial keel



Incision



Anatomical Landmarks



The skin incision should be approximately 1cm medial to the patella. Continue with an arched incision distally to the medial boundary of the patella tendon attachment. Following the skin incision, a medial capsular incision is made (as indicated by the red dotted lines). Expose the patella border and resect osteophytes. To improve visualization, resect the anterior menisci and the infrapatellar fat body until the intercondylar eminence is exposed.

- a) Before starting to prepare for the tibial resection extend the leg fully and mark on the femur with a pen the most anterior wear point between the femur and tibia. This point will be referred to as the tide mark, and represents the anterior boundary of the femoral component.
- b) Also mark the antero medial (or antero lateral) contact area of the femur or tibia in order to assist in defining the position of the tibial sagittal cut.
- c) Prior to the start of the procedure the center of the femoral head may be identified by placing an EKG lead over the femoral pulse just below the inguinal ligament to estimate the center of rotation of the hip. This lead can then be palpated through the drapes and used to assess alignment.

1. Tibial Instrument Assembly



a) Assemble the Ankle Clamp, Tibial Alignment Rod, appropriate Tibial Alignment Guide, and Stylus as shown.

2. Tibial Alignment



- a) With the knee flexed, place the Ankle Clamp around the distal tibia just above the malleoli. Place the stylus in the lowest point on the affected side of the tibial plateau.
- b) The Tibial Alignment Rod should be in line with the anatomic axis in the frontal plane and parallel to the anatomic axis in the sagittal plane. A five degree posterior slope is built into the Tibial Alignment Guide.
- c) An optional 1/8" (3.2mm) pin may be placed percutaneously into the anterior tibia through the Alignment Rod just below the locking knob to help stabilize the extramedullary guide while fine adjustments are made in its alignment.



8000-1040
E/M Tibial Ankle Clamp



8636-0040
Tibial E/M Alignment Rod



Tibial Alignment Guide
8636-0063L LM/RL
8636-0063R RM/LL



8636-0030
Tibial Stylus



1/8" Headless Pin

7650-1038 3.5" Long (4 pack)
7650-1039 2.5" Long (4 pack)




7650-1035
Headless Pin Driver

3. Tibial Alignment Guide






- a) Use 1/8" (3.2mm) pins to secure the Tibial Alignment Guide to the proximal tibia through the lower set of straight holes and through the angled hole furthest from the extensor mechanism.
- b) Remove the Ankle Clamp, Tibial Alignment Rod and Tibial Stylus.


4. Tibial Resection Level



- a) The goal of the tibial resection is to create a space which will accept the thinnest component necessary.
- b) The Tibial Cutting Guides are designed to create a 2, 4, 6 or 8mm resection from the tip of the stylus. Slide the desired thickness Tibial Cutting Guide into the Tibial Alignment Guide. Secure the Cutting Guide by sliding the Retaining Pin through the triangular bracket on the anterior aspect of the Cutting Guide and into the slot in the Tibial Alignment Guide.



5. Tibial Sagittal Cut



- a) Slide the Sagittal Cutting Guide into the Tibial Cutting Guide. Position the Guide to create a resection near the attachment of the anterior cruciate ligament. Lock the Guide into place and use a reciprocating saw against the rounded surface to make the sagittal wall cut to the level of the Tibial Cutting Guide.

Note: Be cautious when making the sagittal wall cut and use retractors to protect surrounding soft tissue and ligamentous structures.

6. Tibial Transverse Cut



- a) Remove the Sagittal Cutting Guide.
 b) Use an oscillating saw to make the transverse tibial resection while being careful not to cut beyond the sagittal wall. The reciprocating saw blade may be left in the cut to assist as a stop for the transverse cut.
 c) Remove the Tibial Cutting Guide, leaving the Tibial Alignment Guide in place.



8636-0096
Sagittal Cutting Guide

Tibial Cutting Guide

8636-0082 2mm
8636-0084 4mm
8636-0086 6mm
8636-0088 8mm



Tibial Alignment Guide

8636-0063L LM/RL
8636-0063R RM/LL

Recommended Equipment

277-96-300 (Stryker®)

Reciprocating Saw Blade
 3.06" long x 0.025" thick

2108-103 (Stryker®)


Sagittal Saw Blade
 3.53" long x 0.035" thick

7. Joint Space Assessment



- a) With the knee flexed 90° insert the appropriate Tensing Positioner into the joint space. Tensing Positioners are available in 4 sizes – 8, 9, 10 and 12mm, which correspond to tibial component thickness. Take the knee through a full range of motion with the Tensing Positioner in place to assess:
- i) Flexion and extension gaps
 - ii) Estimated amount of correction
 - iii) Soft tissue balance.
- b) Once the appropriate Tensing Positioner is identified, remove it from the joint space and set it aside.

8. Femoral Sizing



- a) Based on the pre-operative assessment, select the appropriate Femoral Cutting Guide and assemble it to the Femoral Sizing Handle.
- b) With the knee flexed to 90°, insert the assembly into the joint space. The tongue of the Sizing Handle should be against the posterior condyle and the contoured surface of the Femoral Cutting Guide should be against the distal femur. The proper size femoral component is determined when:
- i) The Femoral Sizing Handle is parallel to the femoral shaft; and
 - ii) The anterior aspect of the Femoral Cutting Guide is at the tidemark; and
 - iii) Sufficient medio-lateral coverage is achieved.



Femoral Tensing Positioner

8636-0076	8mm
8636-0077	9mm
8636-0078	10mm
8636-0079	12mm

8636-0080

Tensing Positioner Handle



8636-0064

Femoral Sizing Handle



Femoral Cutting Guide

8636-0149L	XSmall	LM/RL
8636-0149R	XSmall	RM/LL
8636-0150L	Small	LM/RL
8636-0150R	Small	RM/LL
8636-0151L	Med	LM/RL
8636-0151R	Med	RM/LL
8636-0152L	Large	LM/RL
8636-0152R	Large	RM/LL
8636-0153L	XLge	LM/RL
8636-0153R	XLge	RM/LL

9. Femoral Alignment



- a) Slide the selected Femoral Cutting Guide into the Tensing Positioner previously chosen. Assemble the Distal Femoral Alignment Guide to the Femoral Cutting Guide as shown. Place the assembly into the joint space and pass a long Alignment Rod through the Distal Femoral Alignment Guide. The Femoral Cutting Guide is properly aligned when the following conditions are met:
- i) The knee is flexed to 90° and the anterior aspect of the Femoral Cutting Guide is at the tidemark. This establishes proper flexion/extension position.
 - ii) The Alignment Rod references the femoral head, and is parallel to the femoral axis. This provides alignment with the mechanical axis and varus/valgus position.

10. Distal Alignment Guide



- iii) Sufficient medio-lateral coverage is achieved. The Femoral Cutting Guide profile matches the corresponding implant.
- b) By referencing the flat tibial cut while tensing the joint, the Tensing Positioner provides parallelism between the transverse tibial and posterior femoral resections. The Positioner also sets the A-P position of the component relative to the posterior condyle to create a 6mm posterior femoral resection.
- c) Once properly aligned, pin the Cutting Guide. Note that the two X-holes on the inferior aspect of the Cutting Guide angle downward.



8636-0042

Distal Femoral Alignment Guide



3180-2000

External Alignment Rod
(2 pack)

11. Posterior Femoral Resection and Peg Preparation



- a) Remove the Femoral Tensing Positioner and the Distal Femoral Alignment Guide. Advance the 1/8" Femoral Drill through the smaller drill hole on the fin until it stops.
- b) Prepare for the femoral peg using the large Femoral Peg Drill with Stop.
- c) Resect the posterior condyle by cutting along the bottom surface of the Femoral Cutting Guide, taking care to protect surrounding soft tissue. The resulting femoral resection is 6mm for all sizes.

Note: The femoral peg drill may be left in the cutting guide (as shown) for added stability while making the posterior resection.

12. Fin Preparation



- a) Use an oscillating saw to prepare the fin slot so that it matches the geometry of the component. This is facilitated by inspecting the sagittal view of the femoral trial to estimate the proper depth of the fin slot preparation.

Note: The fin runs the entire length of the component and varies in depth.



8636-0026
1/8" Femoral Drill with Stop



8636-0024
Femoral Peg Drill with Stop


Recommended Equipment
2108-103 (Stryker®)
Sagittal Saw Blade
3.53" long x 0.035" thick



EIUS™ Femoral Trials

6636-3-001	XSmall	LM/RL
6636-3-002	Small	LM/RL
6636-3-003	Medium	LM/RL
6636-3-004	Large	LM/RL
6636-3-005	XLarge	LM/RL
6636-3-011	XSmall	RM/LL
6636-3-012	Small	RM/LL
6636-3-013	Medium	RM/LL
6636-3-014	Large	RM/LL
6636-3-015	XLarge	RM/LL

13. Flexion/Extension Gap Check



- a) An optional flexion/extension gap check may be performed using the Flexion/Extension Alignment Blocks which represent the tibial implant in extension and the sum of the tibial and femoral component thicknesses in flexion.
- b) To assess the tension and alignment of the knee in extension, bring the leg into extension and insert the "extension" side of the Block into the joint space. Long alignment rods may be placed through the Block in both the superior and inferior directions to assess mechanical alignment. Tension and alignment of the knee in flexion may also be examined by taking the knee to 90° and inserting the "flexion" side of the block.

14. Bur Template Placement



- a) Remove the Tibial Alignment Guide.
- b) Select the proper Femoral Bur Template. Place the Template onto the distal femur, referencing the posterior resection, fin slot and 1/8" peg hole.
- c) Using the Femoral Component Impactor, gently tap the template into place. The template is in the proper position when the posterior plate is flush against the posterior resection. Optional Pin holes on the Femoral Bur Template can be used for additional fixation.

3180-2000

External Alignment Rod
(2 pack)



Flex/Ext/Alignment Block

8636-0031	8mm
8636-0032	9mm
8636-0033	10mm
8636-0034	12mm



Femoral Bur Template

8636-0203L	XSmall	LM/RL
8636-0203R	XSmall	RM/LL
8636-0204L	Small	LM/RL
8636-0204R	Small	RM/LL
8636-0205L	Medium	LM/RL
8636-0205R	Medium	RM/LL
8636-0206L	Large	LM/RL
8636-0206R	Large	RM/LL
8636-0207L	XLarge	LM/RL
8636-0207R	XLarge	RM/LL



3140-0000

Eius™ Femoral
Component Impactor

Recommended Equipment

5120-20-347

Stryker® TPS U-Drill fluted
head drum bur Medium (Also
fits Midas Rex and Zimmer
UltraPower)

15. Distal Femoral Burring



- a) Operate the EIUS™ system bur within the guide channel of the Femoral Bur Template. This prepares the distal femur to accept the profile of the femoral component. Proper depth is achieved when the bur stop is flush against the polished surface of the bur template.
- b) Remove the femoral bur template from the femur. Complete the preparation by removing the middle strip of bone using a rongeur or bur.

16. Femoral Trialing and Final Fit



- a) Place the trial component on the femur. Using the Femoral Component Impactor, tap the Femoral Trial until the trial sits fully into the prepared distal femur.



6636-3-XXX

Femoral Trial


3140-0000

Femoral Component Impactor

8636-0070

Femoral Trial Extractor Pliers


17. Tibial Sizing



- Place the Tibial Sizing Guide along the sagittal cut, hooking it over the posterior edge of the tibial resection, note the size indicated.
- Connect the appropriate size keel-less Tibial Trial to the Quick Connect Handle. Place the Tibial Trial on the tibial resection and assess for correct fit. If the margins of the trial extend past the cortical rim, the sagittal cut can be advanced in order to eliminate the overhang.

Note: Size interchangeability between the femoral and tibial components is limited to one up/one down.

18. Trial Reduction



- With the femoral trial and tibial trial in place, take the knee through full ROM. In this way, a check is made to ensure that the components are well centered and that there is no component displacement. The tibial component should be stable and should not lift off or move in the sagittal plane during ROM testing.
- Care should be taken to ensure that there is slight under correction of the overall alignment. Ligamentous tension should be well balanced.



8000-0000

Quick Connect Handle 1/4"



8636-0055


Tibial Sizing Guide




Tibial Trial

6636-3-308	XSmall	8mm	6636-3-510	Med	10mm
6636-3-309	XSmall	9mm	6636-3-512	Med	12mm
6636-3-310	XSmall	10mm	6636-3-608	Large	8mm
6636-3-312	XSmall	12mm	6636-3-609	Large	9mm
6636-3-408	Small	8mm	6636-3-610	Large	10mm
6636-3-409	Small	9mm	6636-3-612	Large	12mm
6636-3-410	Small	10mm	6636-3-708	XLarge	8mm
6636-3-412	Small	12mm	6636-3-709	XLarge	9mm
6636-3-508	Med	8mm	6636-3-710	XLarge	10mm
6636-3-509	Med	9mm	6636-3-712	XLarge	12mm

19. Tibial Keel Preparation



20. Tibial Keel Preparation (continued)














- a) Attach the Alignment Handle to the appropriate Tibial Punch Guide and position on the prepared tibia.
- b) Use the Tibial Impactor on the top rounded edge of the Punch Guide to advance the spikes into the tibia and seat the Guide flush on the tibia. If additional fixation is required, 1/8" headless pins may be drilled into the tibia through the angled anterior holes just above the Alignment Handle.

- c) Place the appropriate Tibial Drill Bushing into the slot of the Punch Guide. Pass the Tibial Keel Drill with Stop through each hole of the Drill Bushing.
- d) Remove the Drill Bushing and impact the appropriate Tibial Punch into the same slot in the Punch Guide until it hits the stop.

Note: A long Alignment Rod can be passed through the Alignment Handle distally as an additional tibial placement check.

	Tibial Punch
	8636-0060 XS/S
	8636-0061 M/L
	8636-0062 XL
	8636-0081
	Eius™ Tibial Impactor
	8636-0025
	Tibial Keel Drill with Stop

	Tibial Punch Guide
	8636-0015L XS LM/RL
	8636-0015R XS RM/LL
	8636-0016L S LM/RL
	8636-0016R S RM/LL
	8636-0017L M LM/RL
	8636-0017R M RM/LL
	8636-0018L L LM/RL
	8636-0018R L RM/LL
	8636-0019L XL LM/RL
	8636-0019R XL RM/LL

	Tibial Drill Bushing
	8636-0020 XS/S
	8636-0021 M/L
	8636-0022 XL
	3180-1000
	External Alignment Handle

21. Final Implantation




- a) Placement of the final implants is facilitated by cementing the tibial component first.
- b) When cementing each component, apply cement to the keel and peg preparations as well as the components. Impactors are provided for both the tibia and femur.
- c) A curette or bent nerve hook may be used to clear out any excess cement paying particular attention to the posterior region of the implant and margins of the tibial eminence.

Note: If sclerotic bone is present, a drill or bur may be used to perforate the bone to improve cement penetration.

Note: Placing a sterile gauze or similar cloth in the posterior joint capsule prior to implantation and slowly removing the cloth after implantation may drag excess cement out of the joint. A dental mirror may also be useful in cement removal.


PREPARATION AND TRIAL TRAYS

TIBIAL PREPARATION TRAY



8636-9001


FEMORAL PREPARATION TRAY



8636-9002


PREPARATION AND TRIAL TRAYS

FEMORAL AND TIBIAL TRIAL TRAY



8636-9003

TIBIAL KEEL PREPARATION TRAY



8636-9004

EIUS Unicompartmental Knee System

IMPLANTS	FEMORAL COMPONENT	TRIALS
LEFT MEDIAL/RIGHT LATERAL	SIZE	LEFT MEDIAL/RIGHT LATERAL
6636-2-001	Extra Small	6636-3-001
6636-2-002	Small	6636-3-002
6636-2-003	Medium	6636-3-003
6636-2-004	Large	6636-3-004
6636-2-005	Extra Large	6636-3-005
RIGHT MEDIAL/LEFT LATERAL	SIZE	RIGHT MEDIAL/LEFT LATERAL
6636-2-011	Extra Small	6636-3-011
6636-2-012	Small	6636-3-012
6636-2-013	Medium	6636-3-013
6636-2-014	Large	6636-3-014
6636-2-015	Extra Large	6636-3-015

EIUS Unicompartmental Knee System

IMPLANTS	TIBIAL COMPONENT	TRIALS
LEFT MEDIAL/RIGHT LATERAL	SIZE	LEFT MEDIAL/RIGHT LATERAL
6636-2-308	Extra Small 8mm	6636-3-308
6636-2-309	Extra Small 9mm	6636-3-309
6636-2-310	Extra Small 10mm	6636-3-310
6636-2-312	Extra Small 12mm	6636-3-312
6636-2-408	Small 8mm	6636-3-408
6636-2-409	Small 9mm	6636-3-409
6636-2-410	Small 10mm	6636-3-410
6636-2-412	Small 12mm	6636-3-412
6636-2-508	Medium 8mm	6636-3-508
6636-2-509	Medium 9mm	6636-3-509
6636-2-510	Medium 10mm	6636-3-510
6636-2-512	Medium 12mm	6636-3-512
6636-2-608	Large 8mm	6636-3-608
6636-2-609	Large 9mm	6636-3-609
6636-2-610	Large 10mm	6636-3-610
6636-2-612	Large 12mm	6636-3-612
6636-2-708	Extra Large 8mm	6636-3-708
6636-2-709	Extra Large 9mm	6636-3-709
6636-2-710	Extra Large 10mm	6636-3-710
6636-2-712	Extra Large 12mm	6636-3-712
RIGHT MEDIAL/LEFT LATERAL	SIZE	RIGHT MEDIAL/LEFT LATERAL
6636-2-318	Extra Small 8mm	6636-3-308
6636-2-319	Extra Small 9mm	6636-3-309
6636-2-320	Extra Small 10mm	6636-3-310
6636-2-322	Extra Small 12mm	6636-3-312
6636-2-418	Small 8mm	6636-3-408
6636-2-419	Small 9mm	6636-3-409
6636-2-420	Small 10mm	6636-3-410
6636-2-422	Small 12mm	6636-3-412
6636-2-518	Medium 8mm	6636-3-508
6636-2-519	Medium 9mm	6636-3-509
6636-2-520	Medium 10mm	6636-3-510
6636-2-522	Medium 12mm	6636-3-512
6636-2-618	Large 8mm	6636-3-608
6636-2-619	Large 9mm	6636-3-609
6636-2-620	Large 10mm	6636-3-610
6636-2-622	Large 12mm	6636-3-612
6636-2-718	Extra Large 8mm	6636-3-708
6636-2-719	Extra Large 9mm	6636-3-709
6636-2-720	Extra Large 10mm	6636-3-710
6636-2-722	Extra Large 12mm	6636-3-712



For cemented use only in the U.S.

Stryker® is a registered trademark of Stryker Corporation.
Eius™ is a trademark of Howmedica Osteonics Corp.

©Howmedica Osteonics 2002 Lit. No. **66364018LI-1** 3M 4/02 MC/ITP Printed in U.S.A.

stryker[®]
Howmedica
OSTEONICS

59 Route 17 South
Allendale, NJ 07401
howost.com